

Tungsten Carbide Cutters | **UM**



High material reduction and smooth surfaces on metal alloys with only one cutter.

The special Komet UM cutting unites three different types of toothing on just one instrument and therefore meets requirements which up to now appeared to be contradictory. The UM cutters offer a range of advantages over tungsten carbide instruments with conventional toothing:

Controlled working by varying the contact pressure

Using UM tungsten carbide cutters, the technician himself decides whether he wishes to place more emphasis on higher material reduction or on smoother surfaces.

High contact pressure = higher material reduction

Low contact pressure = better surface quality

Laboratory tests confirm:

UM cutters produce a better surface than cutters with conventional staggered toothing (E-toothing) and the surface quality achieved is equal to that produced with fine (EF) toothing.

UM cutter				
wwwwwww	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	m	mmm	.
Fine staggered toothing (EF)				
spr	d marine	www.h	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	, m
Staggered tooth	ning (E)			
Mummum	~~~^^	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	mannt	nn,
Comparative dia	agram of th	e surface qu	ality	

(Roughness µm)

Smooth operation with reduced vibration

The completely new blade configuration also contributes to a smoother operation without strain to the operator's wrist.

Economic use

The use of a particularly fine, hot-isostatically pressed tungsten carbide guarantees sharp, unmarred blades and a long service life.

Unique design

Due to its red ring and golden shank, the UM cutter is clearly distinguishable from cutters with conventional toothing.

Application:

1. Reworking of the border line of the subgingival bar/retention lattice with the H77UM.104.023

2. Trimming/smoothing the subgingival bar with instrument H139UM.104.023

3. Shaping of the surface of the mandibular model cast construction, lingual view, with instrument H79UM.104.040

4. Precise working and definition of the margins of a mandibular model cast construction using the H251UM.104.060









Recommendations for use:

- To be used in the laboratory handpiece High contact pressure = higher material reduction Low contact pressure = better surface quality
- Recommended speed: Precious metal O_{opt} 25.000 rpm Non-precious metal and model cast O_{opt} 15.000 rpm



H250UM.104.040