Chirurgia *Surgery*

Rada naukowa:

dr Ivo Agabiti, Pesaro, Włochy
www.sonosurgery.it

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Microsoft Tag Readers

Thanks to the Microsoft Tag Technology, anyone owning a smart phone (Android, Blackberry or iPhone) can play videos featuring selected Komet products.

The service is easy to use: All you need to do is go to your phone's Internet Browser, type in the link <http://gettag.mobi> and follow the directions. Once the software has been uploaded, you can scan the tags.

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Some of the products and designations mentioned in the text are trademarked, patented or copyrighted.

The absence of a special reference or the sign © should not be interpreted as the absence of legal protection.

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This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

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Colours and products subject to alterations. Printing errors excepted.

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4

Extremely fine incisions in bones

Developed in close cooperation with Dr. Ivo Agabiti from Pesaro, Italy, the Sonosurgery sonic tips are a great leap forward in the quality of oral surgery. The innovative tips allow particularly gentle, conservative work on bones within the scope of restorations involving implants.

The Sonosurgery sonic tips SFS100, SFS101 and SFS 102 are indicated for slitting of the crest and lateral incisions in the bone as part of bone spreading treatments and for gentle extraction of a tooth from its alveolar compartment.

The saucer-shaped sonic tips SFS103 and SFS104 as well as the oval SFS105 are particularly suited for gently detaching the sinus membrane within the course of an external sinus lift operation.

For further useful information please visit the website www.sonosurgery.it.

SFS 100/101/102:

- The extremely fine cuts of merely 0.25 mm guarantee maximum conservation of substance in the bone structure.
- The blade length of 10.7 mm provides sufficient penetration depth to allow vertical counter-sinking into the spongy bone; the width of the blade is 3.5 mm.
- Thanks to their oscillating movements, the tips are particularly gentle on the soft tissue.
- The sonic tips are used in oscillating mode in a suitable sonic hand piece, e.g. Komet SF1LM, which is placed on a dental turbine. A special motor is not required. The operating frequency of the oscillating, elliptic movements is approx. 6,000 Hz.
- Sterile cooling solution is supplied via a cooling adapter equipped with small tubes that are screwed in between the sonic hand piece and the sonic tip.

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
Stal nierdzewna

*Sagittal**Cutting thickness 0.25 mm, cutting depth 10.7 mm**Stainless steel***SFS 101**

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Stal nierdzewna

*Axial**Cutting thickness 0.25 mm, cutting depth 10.7 mm**Stainless steel***SFS 102**

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Stal nierdzewna

*Straight**Cutting thickness 0.25 mm, cutting depth 10.7 mm**Stainless steel*

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The sonic tips SFS120, SFS121 and SFS122 facilitate minimally invasive surgical crown extension without the creation of a flap.

These sonic tips were developed in close cooperation with Dr. Thomas Schwenk and Dr. Marcus Striegel, Nuremberg, Germany. They carry out a surgical crown extension for example in cases where the biological width has been violated or when an asymmetric dental arcade or a "gummy smile" are to be rectified.

- *The tip is available in two diameters, allowing an optimal aesthetic treatment of the eye-catching areas. SFS120: Size 020 for teeth 2 and 4, size 030 for teeth 1 and 3.*
- *Posterior use: Even in the posterior region, the bone is often too close to the crown margin.*
- *The treatment is carried out in the tried and tested clinical manner: Having measured the patient's individual biological width, the treatment is planned with the help of a mock-up, wax-up or template. After the subsequent gingival correction, the osteotomy to recreate the biological width without creation of a flap is carried out with the sonic tip by slightly pivoting the tip in mesial and distal direction.*
- *To be used in a suitable sonic hand piece, e.g. Komet SF1LM*

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Stal nierdzewna

Minimally invasive surgical crown extension

Stainless steel

MaxilloPrep Bone

MaxilloPrep Bone

dr Stefan Neumeyer

Bone grafting and congruent augmentation

The MaxilloPrep Bone centring burs TRC and trepan burs TRT are indicated for the safe preparation of bone cylinders. The inventor, Dr. Stefan Neumeyer from Eschlkam, Germany, also uses a combination of these instruments to carry out a systematic and extensive bone augmentation, because the three diameters 045, 055 and 065 are ideally coordinated, allowing the fixation of bone cylinders of matching shapes. Dr. Neumeyer's basic idea is an integral evaluation of both the hard bone substance and soft tissue, in order to guarantee an optimum aesthetical clinical result. The line of MaxilloPrep instruments offers optimally adapted instruments and accessories to suit this comprehensive range of treatment.

- The TRC centring bur creates a circular guide groove. The centring tip of the bur enables exact and safe positioning on the bone substance.
- The TRT trepan bur featuring a state-of-the-art blade configuration, arranged in pairs, is used for the preparation of the bone cylinder. The large perforations in the working part allow both unobstructed view of the operative site and easy removal of the bone cylinders.

The innovative bone chip extractor 9126 as suggested by Dr. Stefan Neumeyer is used for precise, yet gentle removal of bone chips.

- After exposing the bone substance by means of a very thin cut, the bone chip extractor is axially positioned on the bone. The self-centring tip penetrates the bone substance and the funnel-shaped bone chip extractor fills with moist bone chips.

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*Preparation of bone cylinders,
creating a circular guide groove*

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*Creation of the bone bed using centring
burs with increasing or decreasing diameter*

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
Minimally invasive vertical augmentation

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*Removal of bone chips from
spongy bone substance*

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

MaxilloPrep Centring bur for creation of a guide groove caused by centring tip, stainless steel

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

MaxilloPrep trepan bur for the preparation of bone cylinders (transplants), stainless steel

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MaxilloPrep Bonefix

dr Martin Dürholt

Osteosynthesis screws for fixing bone grafts

Prior to carrying out a restoration with implants, the implantologist is frequently faced with a situation where there is not enough bone available. Developed in cooperation with Dr. Martin Dürholt, the MaxilloPrep Bonefix osteosynthesis screws allow exact placing and fixing of the bone graft at the recipient site.

The system comprises high-grade screws made of pure titanium in size 1.5 mm. They are available in lengths of 6, 9 and 12 mm. These osteosynthesis screws are provided with a continuous fast-cutting thread which ends directly at the flat screw head. Secure hold in the handle adapter is guaranteed by the hexagonal reception with a standard width across flats of 1.2 mm. The attachable handle is particularly light in weight and perfectly ergonomic.

Excellent vision during re-entry thanks to the blue anodised surface of the screw.

Advantages

- Small number of components
- "Minitray" suitable for sterilisation
- Protected removal of screws without contamination
- Optimum recess module and high-quality sterilisation container
- Can be used in combination with existing burs
- Compatible with conventional implant screw drivers
- Suitable for all established bone augmentation techniques

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
MaxilloPrep Bonefix

Osteosynthesis screw, hexagonal reception 1.2 mm, pure titanium, blue anodized

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Set sterilisation container A8 (90 x 90 x 55 mm) and insert tray with recesses (without instruments)

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dr Stefan Neumeyer

Bone spreading and bone condensing

Developed in cooperation with Dr. Stefan Neumeyer, the MaxilloPrep Spread and Condense screws made of pure titanium guarantee effective, yet conservative bone spreading and condensing as part of treatments involving implants. These screws mobilize the bone during expansion and facilitate the creation of the bone bed by making use of the residual elasticity of the bone.

- *The Spread and Condense screws are available in 5 sizes. They are screwed in and out in ascending order, until the required target size is reached. The sizes 4 and 5 also come with a triangular cross-section for flexible spreading of extremely hard bones.*
- *Produced according to up-to-date clinical know-how and the latest production techniques, these screws are particularly comfortable to use. Thanks to their double thread, these screws can be screwed in and out in no time at all. Asymmetrical thread flanks with a smaller angle in inserting direction facilitate the insertion of the screws. Penetration into the bone may also be assisted by percussion, i.e. by gently hammering the screw into the bone. To further facilitate the work of the operator, an optional mechanical wrench and a ratchet for use with the manual wrench are also available.*
- *The screws come without o-rings, which makes them particularly easy to clean. A clearly audible clicking sound confirms that the screw has securely locked into the hand wrench.*

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

Bone Spreading with round and triangular screws

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
Stal nierdzewna

MaxilloPrep Spread-Condense

Triangular bur for centering bones, stainless steel

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
Pilot bur for implantology, stainless steel

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Wielkość 5: 4,0-4,8 mm

Czysty tytan

MaxilloPrep Spread-Condense screw

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
4596.000

Set MaxilloPrep Spread-Condense
 Set MaxilloPrep Spread-Condense

558L15.000.1	1	
558L15.000.2	1	
558L15.000.3	1	
558L15.000.4	1	
558L15.000.5	1	
559L15.000.4	1	
559L15.000.5	1	
186A.204.1	1	
210L16.205.013	1	
210L16.205.018	1	
538.000.	1	
538.204.	1	
534.000.	1	
537.000.	1	

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1 Spread-Condense screw, instrument set, sterilisation container and insert tray each

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Angle Modulation System

dr Ernst Fuchs-Schaller

Horizontal widening of the alveolar ridge thanks to flexible mobilization from the inside

The Angle Modulation System developed in close cooperation with Dr. Ernst Fuchs-Schaller allows a minimally invasive horizontal widening of the alveolar ridge even in those cases where there is not enough bone available to successfully insert an implant. The dental surgeon is frequently faced with an initial knife-ridge situation, especially in the lower jaw, in which case the widening of the alveolar ridge is an indispensable precondition.

The new surgical angle modulation method is based on the natural formation of new bone. The basic idea: The cortical lamella is flexibly mobilized from the inside and then successively expanded and erected towards the outside. An essential advantage of this non-traumatic method is that the jaw does not have to be folded up, i.e. the creation of a flap can be avoided. Consequently, the periosteum and the mucosa do not have to be detached and the "biologically active container" remains intact - to promote better healing without complications.

The essential advantage: Expansion of the alveolar ridge and insertion of the implant in just one session.

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4604.000

Set Angle Modulation System A4
Set Angle Modulation System A4

M110.000.016	1		454.000.	1	
M110.000.028	1		186A.204.1	1	
M110.000.034	1		210L16.205.018	1	
M110F.000.028	1		538.000.	1	
M000FC.000.028	1		538.204.	1	
574.000.	1		558L15.000.1	1	
H254E.314.012	1		558L15.000.2	1	
859.314.010	1		558L15.000.3	1	
859.204.010	1		558L15.000.4	1	
210L16.205.013	2		558L15.000.5	1	
			559L15.000.4	1	
			559L15.000.5	1	
			TRC.205.055	1	
			TRT.205.055	1	
			9126.204.042	1	

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
4606.000

Set Angle Modulation System A6
Set Angle Modulation System A6

M110.000.016		1
M110.000.028		1
M110.000.034		1
M110F.000.028		1
M000FC.000.028		1

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

Surgical bone cutters

1. Coarse-grain tungsten carbide:
Porous ingredients; non-calibrated tungsten carbide submitted to a sintering procedure at 1.785 kbar.

2. Blade with nicks:
- produces inaccurate cuts

3. Instrument of average quality:
- traumatic preparation
- inaccurate results
- premature wear
- not cost-effective

4. Fine-grain tungsten carbide:
Sintered according to a HIP process (Hot isostatic pressing), tungsten carbide powder of 1 μ , submitted to a sintering procedure at 2.7 kbar

5. Ideal condition, blades without nicks, gentle, precise cuts

6. Komet quality instrument
- exact preparation
- excellent result
- long service life
- cost-effective work

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Tungsten carbide, ceramics, steel, diamond

Every day, dentists and oral surgeons benefit from the continuous development of enhanced blade configurations and state-of-the-art production methods.

A chronological review takes us back to toothings with continuous blades (H161), followed by the so-called saw teeth, i.e. interrupted blades - the cross-cut (H162) - and finally the "right/right" staggered toothings (H162A) developed by Komet with its extremely high performance and smooth, almost vibration-free operation.

Nowadays, the operator can choose from an extensive range of effective instruments in all shapes and sizes for accurate preparations, cost-effective work and operations with excellent results.

Always pleased to meet new challenges, Komet introduced its range of CeraLine instruments which contains a selection of instruments made of special high-grade ceramics including bone cutters and pilot drills. The raw material used is particularly resistant to exceptionally high loads. This means that these metal-free instruments - that tend to be slightly more brittle due to the nature of the material - can be used in complete safety. The special ceramic used for the white CeraLine instruments is composed of zircon-dioxide ceramic ZrO₂ and aluminium ceramic Al₂O₃ and stabilized by a HIP production process.

A selection of tungsten carbide bone cutters with hard ZrN (Zirconium nitride) coating is available. The ZrN bone cutters are highly resistant to harsh chemicals which might have a damaging effect, for example particularly aggressive chemicals used in instrument baths. The thickness of the coating is merely between 1 and 2 μ m, which is why it does not impair the efficient cutting performance of the instruments.

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

Bone cutter round, with special cutting design for smooth operation

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Bone cutter round, with special cutting design for smooth operation, ZrN coated

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
Combination instrument for conservative preparation of bone tissue and hard tooth substance

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
Bone cutter with high-efficiency cutting design

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Bone cutter, cross cut tapered fissure long

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Bone cutter, round end cross cut tapered fissure

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
Bone cutter, Lindemann, with special staggered tothing

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

Bone cutter, Lindemann, with special staggered tothing, ZrN coated

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Bone cutter, Lindemann short, with special staggered tothing

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



Bone cutter, made of ceramics

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Diamond-coated bone cutters


Compared to cutters, diamond-coated instruments operate in a grinding manner. For example, the round instrument is ideally suited for creating a lateral window during a sinus lift operation.

The Miniflex diamond disc is perfectly adapted for realizing a straight cut to split the bone, even in case of very narrow jaw ridges. To be used with a disc-guard.

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


Miniflex diamond disc for bone-lid method

Apicectomy in the molar area, osteoplastic surgery of the maxillary sinus

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
Miniflex diamond disc for work on bones, ZrN coated

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
Stal nierdzewna

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Pilot bur with internal cooling system (IK) for implantology, stainless steel

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
Stal nierdzewna

Cleaning wire for cleaning the cooling orifices of internally cooled instruments, stainless steel

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Bone cutter with internal cooling system (IK), stainless steel

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Trepan Burs, Stainless steel

The following three trepan burs with different designs are intended for the creation of different kinds of hollow perforations.

227A

Safe explantation

- High-efficiency toothing specially designed for cutting bone, laser depth markings to permit precise exposure of the implants to the required depth, large openings in the working part for a better view of the surgical site.

227B

Preparation of bone cylinders

- High-efficiency toothing specially designed for preparing bone cylinders with a diameter of 2.9 to 8.9, available in different lengths.

227C

Apical ectomy

- Suggested by Dr. Dr. Ayad, Münster, Germany
- Reduced trepanation depth of merely 5 mm.

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

Trepan bur for explantation, stainless steel

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
Trepan bur for preparation of bone cylinders, stainless steel

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Universal pilot burs | Implantology

The universal pilot burs are ideally suited for the initial preparation of the axis and depth of an implant site. If an implant with large diameter is to be placed, it is recommended to enlarge the perforation gradually.

The pilot burs with particularly small diameter can be used for palpating the implant position, perforation of bone cylinders as well as for creating a continuous perforation line during bone spreading.

Advantages:

- Large chip spaces for good chip removal
- Pyramid-shaped, special instrument tip for easy penetration
- Effective cutting
- Lasered depth markings at intervals of 2 mm, starting at 8 mm from the instrument tip
- Size and instrument length lasered on to the shank for easy identification

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





System tulejek do nawierceń

Guide Sleeve System

This guide sleeve system is an easy and cost-effective option for preparing a pilot bore using a template. When used in combination with a template, the system components allow safe paraxial drilling with defined diameter and drilling depth.

The guide sleeve system promotes paraxial drilling without diverging implant axes to guarantee a successful implantation.

- *Drilling of a pilot hole with a defined length of 9 mm and a diameter of 2 mm.*
- *Radiopaque, metallic guide sleeves to be fixed to a drilling template.*
- *Adapter sleeve to be mounted on to the pilot drill, serving as depth stop*
- *Universal use*

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
Stojak na instrumenty

Bur blocks

A range of functional bur blocks is available for cleaning, disinfection, storage and sterilization of dental instruments.

These bur blocks facilitate the cleaning and disinfection of the instruments in the thermo disinfectant, whereas the sterilization takes place in the autoclave. The blocks are provided with silicone plugs to guarantee secure hold of the instruments and to prevent them from falling out.

- Different types of shank can be combined. The bur block fits turbine, hand-piece and contra-angle instruments.
- Equally suitable for organizing preparation sequences, e.g. rotary instruments to perform a sinus floor elevation.
- 9890L4 / Height 40 mm, suitable for example for bone cutters H254E/LE.314 Trepan burs 227A/B/C.204
- 9890L5 / Height 50 mm, suitable for example for bone cutters H141A, H162A, H166A.104 Pilot drills for implantology 210L16/L19/L20.204 and .205
- 9890L7 / Height 70 mm, suitable for instruments with shank 105, long hand-piece

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Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 37 mm

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Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 47 mm

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Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 67 mm

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Bur block made of stainless steel with 7 light blue silicone plugs as universal holders for sonic tips

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Insert trays and Sterilization container

Systematic insert trays and safe sterilization containers

Many systems of the Komet dental range comprise handy insert trays for the reception of instruments and system components. All components are clearly and logically arranged in the correct sequence ready for use during treatment sessions. Both insert tray and instruments are completely integrated in a high-quality sterilization container.

Insert tray

- The insert tray is made of a special high-performance plastic.
- Antimicrobial effect thanks to the continuous release of silver ions
- High purity grade during idle times
- All system components are clearly arranged.

Sterilization container

- Safe sterilization. Does not require sealing in foil
- No seals, no maintenance required. Lasts through more than 2,000 sterilization cycles
- PTFE filter membrane with a particularly small nominal pore size of 0.2 μm - retains even airborne germs and aerosols
- Easy filter change after 150 cycles
- Two or more containers can be connected laterally
- Transparent material - the contents are visible from the outside

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Insert tray MaxilloPrep Spread-Condense, PP with silver ions, antibacterial (without instruments)

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
Sprzęgło

Coupling

Turbine coupling with integrated tube for the external supply of a sterile cooling liquid

The coupling 9969 was specially developed for use with surgical sonic tips. A sterile cooling agent is supplied externally by means of an integrated tube. Thanks to this coupling, a cooling adapter between the sonic hand piece and the sonic tip is no longer required. The water supply of the dental unit is turned off during surgical interventions. The coupling is suitable for sterilisation and can be reprocessed according to a validated method.

With MULTiflex® coupling, suitable for use with the Komet sonic hand piece SF1LM and other air powered air scalers with MULTiflex® coupling. The coupling fits any 4-hole connection.

* MULTiflex® is a registered trademark of the company KaVo Dental GmbH, Biberach

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

Set cooling adapter SF1979 for sonic tips and mounting wrench 566

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Szczotka z metalu
Spare metal brush


9874

Szczotka z nylonu
Spare nylon brush

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Stal nierdzewna

Extension 15 mm for instruments with contra-angle shank

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Stal nierdzewna

Handle for manual use of instruments with contra-angle shank

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Czysty tytan

Hand wheel to be mounted on top the hand wrench, for increasing the torque

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Osteosynthesis screws made of pure titanium

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The bone graft, perforated with the pilot drill 210L16.

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Countersinking of the screw heads - preparation of the perforations with the round bur H141A.104.027.

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The bone graft is screwed into the cortical bone of the recipient site with one of the screws of the MaxilloPrep Bonefix system.

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Minimally invasive widening of the alveolar ridge

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Typical initial situation: Knife-ridge

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Incision and exposure, instrument is slightly tilted in lingual direction, parallel to the defect

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Mobilization of the bone: Deepening the inner side of the buccal lamella and creating vertical incisions at the mesial and distal end of the cortical bone.

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Slow, successive spreading whilst continually checking the flexibility with the flat modulator.

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Erection of the buccal cortical lamella, insertion of the implant during the same treatment session.

Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25 · 32657 Lemgo
Postfach 160 · 32631 Lemgo · Germany

Export:

Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-329
export@kometdental.de
www.kometdental.de

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MULTiflex® is a registered trademark of KaVo Dental GmbH, Biberach