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

Endodontics

Quality Products
Made in Germany

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Some of the products and designations mentioned in the text are trademarked, patented or copyrighted.

The absence of a special reference or the sign ® should not be interpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products subject to alterations. Printing errors excepted.

As at March 2014

2 - 5	环钻术 <i>Trepanation</i>
6 - 8	清除和挖钻 <i>Removal and excavation</i>
9 - 12	扩孔 <i>Reaming</i>
13 - 15	根管口 <i>Canal entrances</i>
16 - 19	SonicLine <i>SonicLine</i>
20 - 27	手动预备 <i>Manual preparation</i>
28 - 52	机械预备 <i>Mechanical preparation</i>
28 - 32	• F360 <i>F360</i>
33 - 40	• EasyShape <i>EasyShape</i>
41 - 47	• AlphaKite <i>AlphaKite</i>
48 - 50	• EndoPilot <i>EndoPilot</i>
51 - 52	• E-Drive <i>E-Drive</i>
53 - 55	Endo Rescue Kit 取根管断针套装 <i>Endo Rescue Kit</i>
56 - 61	保存和循环再处理 <i>Storage and reprocessing</i>
62 - 70	填充 <i>Obturation</i>
71 - 72	配件 <i>Accessories</i>
73	根管桩 <i>Root posts</i>
74	复合树脂系统 <i>Composite System</i>

Trepanation

Trepanation is always the first step of the sequence of an endodontic treatment. The access to deeper regions in the tooth is created by penetrating into the various existing substances: through dental enamel, metal crowns, crowns with ceramic veneers or zirconium oxide crowns. A crown cutter with a powerful cutting performance is used in those cases where the crown has to be renewed.

A range of specially developed instruments for quick, effective removal of all types of material is available. Each material is cut according to its special properties, especially ceramics and bonded materials such as metal crowns with ceramic veneers.

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
H140.314. ...

016

⊙_{max.} 120000 rpm
用于环钻金属全冠的麻花钻
Twist drill for trepanation of metal crowns

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

For porcelain-fused-to-metal crowns and metal crowns
Apply crown remover at an angle of 45°
(For extremely hard ceramics, as for example zirconium,
use 4ZR-Diamond)

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
*For porcelain-fused-to-metal crowns and metal crowns
Apply crown remover at an angle of 45°
(For extremely hard ceramics, as for example zirconium,
use 4ZR-Diamond)*

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

*For metal crowns
Apply crown remover at an angle of 45°*

●● H 34

		5	5
尺寸 · Size	∅ 1/10 mm	010	012
L	mm	2,0	2,0

FG车针 (FG) FG · Friction Grip (FG)

500 314 138008 ...

●● H34.314. ... 010 012

用于金属全冠在45°角度应用破冠车针
For metal crowns
Apply crown remover at an angle of 45°

H 40

		5
尺寸 · Size	∅ 1/10 mm	012
L	mm	4,0

FG车针 (FG) FG · Friction Grip (FG)


500 314 139008 ...

H40.314. ... 012

用于金属全冠在45°角度应用破冠车针
For metal crowns
Apply crown remover at an angle of 45°

●○ ZR 6801

		5	5	5
尺寸 · Size	∅ 1/10 mm	010	014	018

FG车针 (FG) FG · Friction Grip (FG)

●○ ZR6801.314. ... 010 014 018

球形, 用于环钻氧化锆全冠
Round, for trepanation of zirconia crowns

new

○ 4 ZR


		5	5
尺寸 · Size	∅ 1/10 mm	012	014
L	mm	8,0	8,0

FG车针 (FG) FG · Friction Grip (FG)

○ 4ZR.314. ... 012 014

用于切割氧化锆全冠和桥体
For cutting of zirconia crowns and bridges

Removal and excavation

The second step within the course of an endodontic treatment is the removal of carious lesions and fillings. Round burs operating with low vibration are predominantly used for these tasks. To this end, Komet has set a standard by developing high-end tungsten carbide instruments provided with a unique hybrid toothing with a tip transversing blade and staggered toothing on each blade as well as round burs made of high performance ceramics which allow particularly gentle, intuitive excavation. Sharp blades and an absolute minimum of vibration during use enable the dentist to quickly remove hard dental substance whilst offering the patient the greatest possible comfort during the treatment.


Amalgam and composite fillings can be efficiently removed with specially developed instruments.

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
Round bur, made of ceramics
High-efficiency cutting design for excavating
Slim neck for improved vision

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
Reduced vibration H1SE/
H1SEM compared to
conventional round burs

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

Round

Staggered toothing with high-efficiency cutting design for
excavating

窝洞预备
圆钻
Cavity preparation
Round

8

H 1 SM

		5	5	5	5	5	5	
尺寸 · Size	∅ 1/10 mm	010	012	014	016	018	021	023
美国编号		2SM	3SM	4SM	5SM	6SM	7SM	8SM
FG车针 (FG) FG · Friction Grip (FG)								
	H1SM.314. ...	-	-	014	-	-	-	-
RA车针 (RA) RA · Right-angle (RA)								
	H1SM.204. ...	010	012	014	016	018	021	023
长柄RA车针 (RAL) RAL · Right-angle long (RAL)								
	H1SM.205. ...	010	-	014	-	018	-	023

■ = max. 100000 rpm
■ = max. 160000 rpm

圆钻
高效切割设计，用于挖钻细颈用于改善视野
Round
High-efficiency cutting design for excavating
Slim neck for improved vision

H 32

高效工作的银汞清除车针
Amalgam remover
for efficient work

		5
尺寸 · Size	∅ 1/10 mm	012
L	mm	4,2
FG车针 (FG) FG · Friction Grip (FG)		
	H32.314. ...	012

银汞清除车针
Amalgam remover

5985

清除陈旧的填充物
复合树脂清除车针
Removal of old fillings
Composite remover

		5
尺寸 · Size	∅ 1/10 mm	012
L	mm	4,2
FG车针 (FG) FG · Friction Grip (FG)		
	5985.314. ...	012

实用新型专利，专利 / Utility model, patents
DE 199 08 507
EP 1 031 325

复合树脂清除车针
Composite remover

Apical problems solved the coronal way!

According to this motto for a successful endodontic treatment, the roof of the pulp chamber is completely removed. To create a straight access cavity to the root canal, lateral gain of space and the removal of overhanging material are particularly important. Both abrasive and cutting instruments are available for this purpose. They can be used individually or in combination, as required.

383, 389: Combination instruments with a diamond coated working part and a cutting tip made of tungsten carbide, for lateral substance removal, creating lateral space and gaining depth, if required.

Instruments with safe, non-cutting ends: Diamond and tungsten carbide burs with safe ends are available to protect the base of the pulp chamber and the canal entrance.

15802: With these special tapered diamond instruments, the pulp chamber can be reached with ease.

根尖问题破坏了根冠的通道 (Apical problems solved the coronal way)

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
FG车针 (FG) FG · Friction Grip (FG)

H269QGK.314. ... 016

∅_{max.} 160000rpm
顶端安全无刃，带交叉纹路的牙髓腔扩孔钻
Endo reamer with safe end, with cross cut

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
Combination bur with diamond coated working part and
carbide tip, for lateral substance removal when preparing
an access cavity

10

15802

			5
尺寸 · Size	∅ 1/10 mm		014
L	mm		10,0

FG车针 (FG) FG - Friction Grip (FG)

15802.314. ...

014

⊖_{max} 160000 rpm

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
For creation of an endo access cavity

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


Canal entrances

Established Reaming made of stainless steel or tungsten carbide are used for the retrieval of root canals and tapered reaming in the straight, coronal part of the root canal. The endodontic tungsten carbide bur H1SML is particularly suited for the creation of isthmuses. With its help, an isthmus can be created along a part or the entire length to detect and expose even deeply hidden canals.

Round pulp burs help to penetrate the canal entrance. For the subsequent reaming of the straight, coronal portions of the root canal, use Gates-Glidden burs, applying the Crown-Down technique. Alternatively, burs with a large taper are equally suited for this purpose. Thanks to the generously enlarged canal entrance, all subsequent instruments can work in the deeper regions of the root canal with ease.

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
⊙_{max.} 20000 rpm
不锈钢G型扩孔钻
Reamer Gates Glidden "G", stainless steel

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

Reamer Gates Glidden "G", short, stainless steel

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For isthmus preparation, slim neck for improved vision

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


Reamer “B1”, stainless steel

SonicLine

Orthograde preparation of the pulp chamber and preparation of the cervical third of the root canal as part of an endodontic treatment

Advantages:

- Quick preparation and removal of old root fillings
- Easier retrieval of root canals
- Enlargement of obliterated canals
- Preparation of straight canal access cavities without weakening the crown
- Useful for removing hard root fillings, cements or posts
- Controlled, gentle preparation without steps and protruding material
- Excellent vision

The sonic tips can be used in the air driven sonic hand piece SF1LM (air scaler) and reprocessed with the help of a rinse adapter in a Miele washer/disinfector.

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大圆头
初步预备牙髓腔通道并清除突出材料
*Large bud
Initial preparation of the access cavity and removal of protruding substance*

SF 67

		1
L	mm	10,0
角度 · Angle	α	125°
SF67.000. ...		•

锥形头
修复细小和钙化的根管，在修复时打开根管的上部
Tapered
Retrieval of fine and calcified canals, opening of the upper canal portions during revision

SF 68

		1
L	mm	10,0
角度 · Angle	α	112°
SF68.000. ...		•

锥形头
替代带有更大角度的SF67
Tapered
Alternative to the SF67 with a more pronounced angle

SF 69

		1
L	mm	6,0
SF69.000. ...		•

小圆头
精修牙髓腔通道，最小程度复新牙本质层，清除残留的根管填充物
Small bud
Finishing of the access cavity, minimal refreshing of the dentin layer and removal of residues of root canal fillings

SF 70

		1
角度 · Angle	α	122°
L	mm	10,0
SF70.000. ...		•

锥形头
扩张较长和较宽的根管，清除折断的器械，清除由牙胶和软粘剂构成的根管填充物
Tapered
Enlargement of long and wide canals, unblocking of fractured instruments, removal of root fillings made of gutta-percha and soft cements

Sonic tip for activating rinsing liquids during endodontic treatments

Thanks to sound activated movements and micro currents, the efficiency of rinsing liquids inside the root canal is greatly increased, safely removing bacteria, residues of pulp tissue, loose dentin chips and the smear layer.

Advantages:

- Even more efficient rinsing of the root canal
- Made of highly flexible nickel titanium, with a titanium-nitride coating
- Small instrument diameter for rinsing narrow root canals
- The same instrument diameter for all canals, thus eliminating the need to change the instrument
- Safe end and absence of toothing to avoid inadvertent removal of substance from the canal wall
- Laser marks to indicate the depth

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
9952.000

尺寸 · Dimensions mm 90 x 65 x 22

含有7个淡蓝色硅胶塞的不锈钢车针架，作为声波工作尖的通用型支架
Bur block made of stainless steel with 7 light blue silicone plugs as universal holders for sonic tips

Three basic types of files have been used for decades for the manual preparation of root canal systems: K-burs, K-files and Hedstroem (H) files. The preparation of root canals is ruled by the standard ISO 3630 for root canal instruments. In compliance with ISO, the taper of manual files is 2% and all files have a toothed working part with a length of 16 mm.

Beyond the standard determined by ISO, the manual preparation instruments made by Komet offer the following advantages:

- Radiopaque silicone stoppers preassembled at our works
- Special ergonomic handles to protect from horizontal or vertical slipping
- The handle is provided with a hole to insert a safety thread
- Ground depth marks at 18 - 19 - 20 mm from the instrument tip
- NiTi files are provided with a gold coloured titanium nitride coating to protect the file against premature wear caused by the sterilisation procedure

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
9107.634. ...

020 025 030 035 040 050 060

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Nerve broaches, stainless spring steel

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
Instruments for root canal preparation, root canal

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
Instruments for root canal preparation, root canal reamers “K”, stainless spring steel

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
NiTi-Burs “K”, nickel-titanium alloy, TiN-coated

171 FC 21
171 FC 25

Sortimente:
Assortments:

L = 25 mm

175(25).654.S0

2 x 006 - 010

17525

		6	
尺寸 · Size	Ø 1/100 mm	010	
L	mm	25,0	
手柄 · Handle			
	17525.654. ...		010

用于探查和检查开放根管的手用锉，不锈钢，热处理钢
Manual file for probing and checking patency of root canals, stainless, heat-strengthened steel

Sortimente:
Assortments:


L = 21 mm

173(21).654.S1

1 x 015 - 040

173 25.654.S1

17321
17325
17328
17331


		6	6	6	6	6	6	6	6
尺寸 · Size	Ø 1/100 mm	006	008	010	015	020	025	030	035
手柄 · Handle									
	340 654 645452 ...	17321.654. ...	006	008	010	015	020	025	030
	340 654 646452 ...	17325.654. ...	006	008	010	015	020	025	030
	340 654 647452 ...	17328.654. ...	006	008	010	015	020	025	030
	340 654 648452 ...	17331.654. ...	006	008	010	015	020	025	030

用于根管峡部的预备，K型根管锉，不锈钢簧钢
Instruments for root canal preparation, "K" files, stainless spring steel

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
Instruments for root canal preparation, “K” files, stainless

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
NiTi-Files “K”, nickel-titanium alloy, TiN-coated

17421
17425
17428
17431

Sortimente:
Assortments:

L = 21 mm

17421.654.S1

1 x 015 - 040

17425.654.S1


17428.654.S1

尺寸 · Size	Ø 1/100 mm	008	010	015	020	025	030	035	040
手柄 · Handle									
340 654 650453 ...		008	010	015	020	025	030	035	040
340 654 651453 ...		008	010	015	020	025	030	035	040
340 654 652453 ...		008	010	015	020	025	030	035	040
340 654 653453 ...		008	010	015	020	025	030	035	040

根管预备器械, “H”型根管锉, 不锈钢
Instruments for root canal preparation, Hedstroem files
„H”, stainless spring steel

17421
17425
17428
17431

Sortimente:
Assortments:

L = 25 mm


17425.654.S2

1 x 045 - 080

尺寸 · Size	Ø 1/100 mm	045	050	055	060	070	080	090	100	110
手柄 · Handle										
340 654 650453 ...		045	050	055	060	070	080	090	100	110
340 654 651453 ...		045	050	055	060	070	080	090	100	110
340 654 652453 ...		045	050	055	060	070	080	090	100	110
340 654 653453 ...		045	050	055	060	070	080	090	100	110

根管预备器械, “H”型根管锉, 不锈钢
Instruments for root canal preparation, Hedstroem files
„H”, stainless spring steel

174 FC 21
174 FC 25

Sortimente:
Assortments:

L = 21 mm

174FC21.654.S1

1 x 015 - 040

174FC654.S2

1 x 045 - 080

		6	6	6	6	6	6	6	6	6	6	6	6
尺寸 · Size	Ø 1/100 mm	015	020	025	030	035	040	045	050	055	060	070	080


手柄 · Handle

174FC21.654. ...	015	020	025	030	035	040	045	050	055	060	070	080
174FC25.654. ...	015	020	025	030	035	040	045	050	055	060	070	080

NiTi "H"型锉, 镍钛合金, TiN涂层
NiTi-Hedstroem files "H", nickel-titanium alloy, TiN-coated

F360 – Anything other than complicated

- Root canal preparation with only 2 files
- Each file is used in rotary mode to its full working length
- Single-use files
- Sterile-packed

Applied intelligence

Prior to using the F360 files, an Opener (e.g. AK10L19) is used to ream the canal access zone to a larger size and remove all infected substance. Subsequently, patency of the canal has to be ensured by means of suitable manual instruments. In most cases, the root canal can now be prepared with just 2 files. The red F360 file in size 025 is first in line, followed by the green F360 file in size 035. The use of the first F360 file in size 025 may be sufficient in case of very narrow canals. On the other hand, the use of additional files in sizes 045 (white) and 055 (red) may be required in case of wide canals.

The F360 files are used in rotary “picking” motion, i. e. they are moved along the entire working length in a dabbing, picking manner so that the files are kept moving at all times.

Thanks to their taper of 04, all F360 files are extremely flexible so that they are perfectly able to adapt to different root canal morphologies.

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
Opt. 250 - 350 rpm
扭力: 1,8 Ncm

Well thought-out in every detail

- Manageable file sequence: only 2 files are needed to suit most canal anatomies
- Use in rotary mode in a conventional torque-limited endodontic contra-angle or motor
- All files are used with an identical torque of 1.8 Ncm so that there is no need to adjust the torque of the endodontic contra-angle when changing the file
- Simple use according to the single-length principle
- Optimum respect of the original canal anatomy
- Excellent cleaning result and quick mechanical preparation thanks to the innovative instrument design:
 - > cross-section in the shape of a double S for increased flexibility, large chip spaces to flush out all debris, blade twist with dynamic increase.
- Reduced risk of instrument fracture and no risk of cross-contamination as the F360 files are intended for single use only
- No need to clean, disinfect and sterilise instruments. No need to document and control the number of uses as the F360 single-use files are sterile packed.
- Flexible files made of NiTi with small taper 04 assure optimal adaptation to nearly all canal anatomies
- The files are also available in the additional sizes 045 and 055 and come in 3 lengths (L21, L25 and L31) to allow for optimal preparation of wider canals

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
F360 Fill

F360 Fill is a carrier-based filling system for thermoplastic, three-dimensional root canal fillings. The F360 Fill obturators consist of a plastic core coated with thermoplastic gutta-percha. The obturator is heated in the F360 Fill oven in order to guarantee a tight obturation of the root canal. F360 is a universal post system that ideally complements for example the F360 files.

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F360 Fill Verifier

Instruments to determine the correct size of the F360 Fill

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F360 Fill Obturators

Obturator consisting of a plastic core, coated with thermo

plastic gutta-percha, to be heated in the F360 Fill Oven to


achieve a tight, three-dimensional filling

Universal obturator system, suitable for example for F360

files

Color coded and radiopaque

AlphaKite and EasyShape

are two systems of files made of nickel titanium designed for the mechanical preparation of the root canal. AlphaKite is used in the Crown-Down technique, whereas the EasyShape files are used at the full working length. Both systems share the following features:

Instrument sizes – ISO colours

The size of the files can be recognised by a coloured section on its shank. The sizes were allocated according to the ISO colour scheme, according to which a file with a yellow mark corresponds to size 020, a file with a red mark means size 025, blue equals size 030, green 035 and black 040.

Taper – Laser mark on the shank

Apart from the coloured section, the size of the EasyShape files can be identified by a laser mark on their shank. The laser mark shows the figure number and the taper of each file, e.g. ES04.105 white describes the EasyShape file with a taper of 04 in size 015.

TiN plated

All files of the EasyShape- and AlphaKite-Systems are coated with a special titanium nitride layer. Apart from the special blade configuration, this surface coating is the reason for the excellent, long lasting cutting performance of the instruments. What's more, it prevents premature blunting caused by sterilization or contact with sodium hydrochloride.

Controlled work length – ground-in depth marks, preassembled silicone stoppers

In order to stick to the work length of the instruments irrespective of the work conditions, the files have been provided with ground-in depth marking rings (at 18, 19, 20 and 22 mm from the instrument tips) and preassembled, radiopaque silicone stoppers.

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
EasyShape

The most outstanding feature of the EasyShape system is its very sharp and effective blade configuration in the shape of a double-S. The EasyShape system was developed in close collaboration with Prof. Dr. Edgar Schäfer of the University Hospital in Münster, Germany.

Great cutting force

The principle of the EasyShape file system is to work with files along the entire length of the root canal. Each EasyShape file is moved up and down the canal wall in short, continuous movements. Each file is used to its full length, provided that there is no resistance in the canal. The files are used in "picking motion", i.e. they are moved along the canal wall in a dabbing, picking manner. To prevent the particularly sharp toothing of the instruments from jamming in the canal wall, the files are kept moving at all times.

The files are provided with non-cutting tips, which is why the root canals have to be perfectly free of obstructions prior to the start of the preparation. Complete patency of the canal must be guaranteed because in general, EasyShape instruments such as NiTi files are unable to penetrate blocked root canals. Guide the EasyShape files through the canal in the usual manner without excessive force. The generous chip spaces of the files allow effective evacuation of debris from the canal.

- Preparation along the entire working length
- Manageable system containing 6 files
- Double-S toothing for outstanding cutting force
- Titanium nitride coating for enduring cutting performance

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
Manageable system


The EasyShape files are used in ascending order. The straight canal access zone can be reamed to a larger size at the start of the treatment, if required. To this end, the access reaming file AK10L15.035 with a taper of 10 from the KOMET AlphaKite system can be used.

The file ES04.015 is first in line, i.e. a size 015 file with a taper of 04. This is followed by inter-coordinated files 05.020, 06.025, 05.030, 04.035 up to size 04.040. The files are used in ascending order; the system is self-explanatory and consists of only a few steps. The instruments are always used in the same sequence - starting with size 015 until the required apical target size has been reached. In case of narrow canals, this can already be achieved with the fourth file, i.e. taper 05.030 blue. The instruments are available in three different lengths: L21, L25 and L31.

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


EasyShape file taper 04, spiral-shaped colour code
File for canal preparation to the full working length in picking motion, nickel-titanium with TiN coating

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
EasyShape file taper 05, spiral-shaped colour code
File for canal preparation to the full working length in picking motion, nickel-titanium with TiN coating

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
EasyShape file taper 06, spiral-shaped colour code
File for canal preparation to the full working length in picking motion, nickel-titanium with TiN coating

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

AlphaKite File Taper 10, length 15 mm or 19 mm
Universal reamer for the straight root canal entrance area, nickel-titanium with TiN coating

4584.000

4585.000

EasyShape器械基本套装
EasyShape instrument starter set

EasyShape器械基本套装
EasyShape instrument starter set

●	AK10L19.204.035	2	
○	ES04L25.204.015	2	
●	ES05L25.204.020	2	
●	ES06L25.204.025	2	
●	ES05L25.204.030	2	
●	ES04L25.204.035	2	
●	ES04L25.204.040	2	

2x7支根管锉，工作长度25mm
2 x 7 instruments, for a working length of 25 mm

●	AK10L19.204.035	1	
○	ES04L25.204.015	1	
●	ES05L25.204.020	1	
●	ES06L25.204.025	1	
●	ES05L25.204.030	1	
●	ES04L25.204.035	1	
●	ES04L25.204.040	1	

7支根管锉，工作长度25mm
7 instruments, for a working length of 25 mm

4620.000

EasyShape根管锉套装
EasyShape instrument set

	AK10L19.204.035	1		
	ES04L2.1.204.015	1		
	ES05L2.1.204.020	1		
	ES06L2.1.204.025	1		
	ES05L2.1.204.030	1		
	ES04L2.1.204.035	1		
	ES04L2.1.204.040	1		

7支根管锉，工作长度21mm
7 instruments, for a working length of 21 mm

4621.000

EasyShape根管锉套装
EasyShape instrument set

	AK10L19.204.035	1		
	ES04L31.204.015	1		
	ES05L31.204.020	1		
	ES06L31.204.025	1		
	ES05L31.204.030	1		
	ES04L31.204.035	1		
	ES04L31.204.040	1		

7支根管锉，工作长度31mm
7 instruments, for a working length of 31 mm

533.000

含银离子聚丙烯塑料制成的EasyShape插入盘, 抗菌 (无器械)
Insert tray EasyShape, PP with silver ions, antibacterial (without instruments)

532.000

尺寸 · Dimensions mm 90 x 90 x 55

灭菌容器, 未密封, 无需维护, 带有150次灭菌周期的灭菌过滤器, 可叠放, 可连接, 由透明UUPS塑料制成
Sterilisation container, no seals, no maintenance, with sterilisation filter for 150 sterilisation cycles, stackable, connectable, transparent PPSU plastic

4588.000

EasyShape套装灭菌容器和插入盘 (无器械)
Kit EasyShape, sterilisation container and insert tray (without instruments)

594.000

含银离子聚丙烯塑料制成, 适用于各种器械长度的EasyShape插入盘, 抗菌 (无器械)
Insert tray EasyShape for all instrument lengths, PP with silver ions, antibacterial (without instruments)

		灭菌容器 Sterilisation container
532.000.	1	EasyShape插入盘 EasyShape insert tray
533.000.	1	EasyShape插入盘 EasyShape insert tray


4619.000

适用于各种工作长度器械的EasyShape套装，含灭菌容器和插入盘（不含器械）
Set EasyShape for all instrument lengths, sterilization container and insert tray (without instruments)

532.000.	1		灭菌容器 Sterilisation container
594.000.	1		适用于各种工作长度器械的 EasyShape插入盘 Insert tray EasyShape for all instrument lengths

AlphaKite

AlphaKite is a system of NiTi files designed for use according to the Crown-Down technique. The files contained in the system can be combined in many different ways to suit the operator's requirements.

Variable

The operator can choose from 6 different tapers to ensure effective mechanical cleaning of the root canal. All AlphaKite files are designed for a working length of 25 mm. Intended for initial reaming of the straight upper canal access area, the files with a taper of 10 come in lengths of 19 mm and 15 mm because they are not intended for work down to the full depth of the canal. Any Crown-Down sequence can be preceded by the reaming of the canal access area with a file of taper 10.

AK10, length 15 mm and 19 mm: universal reamer for the straight canal access area.

AK 08, AK 06: for the upper third of the canal and the initial part of the curve.

AK04: for the central and the apical region of the canal.

AK03.015: this instrument is designed for probing.

AK02: a file with a particularly slim taper for the apical region, notably in cases of extremely curved canals, within the scope of treatments according to the Crown-Down technique.

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

红色顺序——一般根管
Red sequence - average canals

蓝色顺序——较宽根管
Blue sequence - wide canals

Crown-Down technique
Coronal reaming of the canal entrance

First of all, the straight access area of the root canal is prepared with a file of taper 10. Next, a manual file (size 015) is used to create a smooth path for the subsequent files and to determine the working length. The next step is the preparation of the canal with instruments of the apical target size (depending on the size of the file that created the glide path), beginning with a taper of 06 and followed by a taper of 04 up to approx. 3 mm before the end of the preparation. The canal is now prepared to the full apical depth with a file of taper 02.

The canal can subsequently be enlarged by one size with the help of the next larger file with the same taper 02.

In some cases, depending on the curve of the canal, the file with a taper of 04 can be used all the way down to the full working depth, provided that the canal is free of perceptible obstructions.

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

Safe and effective

All AlphaKite files are provided with the typical Komet kite-shaped cross-section. One active and three supporting cutting angles ensure an effective and safe preparation. AlphaKite files are designed for use in a torque-limited motor. The files glide into the canal. Large chip spaces guarantee optimum evacuation of debris from the root canal. The effective main cutting angle guarantees the removal of the correct amount of substance. Three further blades ensure an optimum interplay of the rotational forces. The AlphaKite files are easy to control.

During practical use, allow the instruments to gently glide into the root canal with intuition. The blade twist is subject to a slight dynamic increase from the tip of the instrument to the end of the working part. This allows any debris to be effectively removed from the canal, thus ensuring an optimum mechanical cleaning result.

- Work according to the Crown-Down technique
- 25 files for versatile use
- Kite-shaped cross-section for safe, effective preparation
- Titanium nitride coating for enduring cutting performance

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

AlphaKite File Taper 10, length 15 mm or 19 mm

Universal reamer for the straight root canal entrance area, nickel-titanium with TiN coating

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
AlphaKite file taper 08, length 25 mm, spiral-shaped colour code

File for initial canal access reaming, nickel-titanium with TiN coating

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

AlphaKite file taper 04, length 25 mm, 2 rings

File for canal preparation, nickel titanium with TiN coating

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
AlphaKite file taper 03, length 25 mm, spiral-shaped colour code

File for initial probing, nickel-titanium with TiN coating

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AlphaKite file taper 02, length 25 mm, 3 rings

File for canal preparation, nickel-titanium with TiN coating

540.000

含银离子聚丙烯塑料制成的AlphaKite插入盘，抗菌（无器械）
Insert tray AlphaKite, PP with silver ions, antibacterial (without instruments)

539.000

尺寸 · Dimensions mm 90 x 90 x 55

AlphaKite灭菌容器，未密封，无需维护，带有150次灭菌周期的灭菌过滤器，可叠放，可连接，由透明PPSU塑料制成
Sterilisation container AlphaKite, no seals, no maintenance, with sterilisation filter for 150 sterilisation cycles, stackable, connectable, transparent PPSU plastic

46


4579.000

AlphaKite 套装灭菌容器和插入盘（无器械）
AlphaKite kit, sterilisation container and insert tray (without instruments)

		☞
539.000.	1	AlphaKite灭菌容器 Sterilisation container AlphaKite
540.000.	1	AlphaKite插入盘 Insert tray AlphaKite

4574.000

4575.000

AlphaKite器械基本套装
AlphaKite instrument starter set


AlphaKite器械套装+定序器
AlphaKite instrument set + sequencer

●	AK10L19.204.035	3	
●	AK06L25.204.020	1	
●	AK04L25.204.020	1	
●	AK02L25.204.020	1	
●	AK02L25.204.025	2	
●	AK06L25.204.025	1	
●	AK04L25.204.025	1	
●	AK02L25.204.030	2	
●	AK06L25.204.030	1	
●	AK04L25.204.030	1	
●	AK02L25.204.035	1	
○	AK03L25.204.015	1	

2x8支根管锉, 工作长度25mm
2 x 8 instruments, for a working length of 25 mm

●	AK10L19.204.035	1	
●	AK06L25.204.020	1	
●	AK04L25.204.020	1	
●	AK02L25.204.020	1	
●	AK06L25.204.025	1	
●	AK04L25.204.025	1	
●	AK02L25.204.025	1	
●	AK02L25.204.030	1	

工作长度为25mm的8支根管锉和9870型定序器
8 instruments, for a working length of 25 mm and sequencer 9870

EndoPilot

The EndoPilot is an endodontic motor, an apex locator and a tool for obturations all in one, ensuring an efficient and safe preparation of the root canal.

This torque and speed controlled endodontic motor is provided with coloured LED lights indicating the direction of rotation, the torque limit or the position of the apex. These features ensure an efficient preparation of the root canal. Thanks to the fully insulated electric contra-angle, the integrated apex locator allows an exact, real-time determination of the length because the operator is always in full control of the actual position of the file.

The EndoPilot is provided with a file library containing the characteristics of the F360, EasyShape and AlphaKite files and other commonly used file systems. The "MyFile" function enables the operator to develop individually adapted sequences and to save them.

A file management tool monitors the wear of each file by means of important wear parameters such as speed, torque and duration of use.

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
EndoPilot
© 410243 | © 410244

To finish off the treatment, a three dimensional, tight filling can be achieved with the warm, vertical condensation technique. The DownPack hand piece is used for filling the apical third, whereas the central and upper third are thermo plastically filled with the BackFill gun.

The EndoPilot is provided with a large, easy to read touch screen which guides the operator through the menu. All functions are found quickly and easily.

The EndoPilot is designed according to a modular construction. The system can therefore be expanded at any time.

Advantages:

- Endodontic motor, apex locator and thermo plastic obturation system all in one
- Fully insulated motor and contra-angle
- File library with many preset file systems
- Option to develop an individual sequence (MyFile)
- Torque limiting function, equally effective in reverse operation
- Apex locator with impulse measurement
- The preparation length is manually variable at the apex locator
- Once the apex has been reached, the motor stops instead of switching to reverse operation, which would transport the chips towards the apex
- Torque reduction in the vicinity of the apex
- Patented file management tool to monitor the wear of the file
- Modular construction of the system
- Clearly laid out menu with easy-to-use touch screen
- Can be updated to allow for future technical developments

An alternative "comfort" version of the EndoPilot is also available. This version comes without the thermo plastic obturation system, however, the obturation module can be fitted at any time.

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
EP0012 EndoPilot comfort plus
EP0011 EndoPilot comfort (没有填充模块)

E-Drive – Torque limited endodontic contra-angle

The E-Drive contra-angle is directly placed on the coupling of the micro motor. It can be used with all commonly used nickel titanium file systems (e.g. F360, Alpha-Kite and EasyShape). Due to the fact that the torque can be transmitted at 5 different levels and the E-Drive can be directly connected to an apex locator, the mechanical preparation of the root canal can take place in complete safety. Should the file jam in the canal because the torque has been exceeded, it can be released by retro rotation (left-right movements).

Advantages:

- The E-Drive can be placed directly onto the coupling of the micro motor
- The torque can be transmitted at 5 different levels
- Optional setting of intermediate levels
- The torque can be set directly at the selection ring of the contra-angle
- The E-Drive can be directly connected to an apex locator
- The apex can be located without using a file clamp
- Transmission 115 : 1
- Maximum speed of the motor: 40,000 rpm
- Performs full rotations 360°
- Safety thanks to retro rotation
- Small head to ensure unobstructed view
- Can be sterilised in the autoclave at up to 135°C

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9938.000

E-Drive
扭力限制的根管治疗弯机
E-Drive
Torque limited endodontic contra-angle

The E-Drive can be directly connected to an apex locator by means of an apex clip, thus eliminating the need for a file clamp. There are two apex clips available for the E-drive hand-piece, allowing the connection of a measuring cable for an endometric measuring device for exact determination of the position of the file in the root canal. The apex locator is attached to the coupling case, the cable of the apex locator is positioned on the clip and the mucosa electrode is attached to the lip.

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The apex clips are available in two different diameters (2.0 mm and 1.5 mm) for connection to various endometric devices.

9939.000

E-Drive根尖夹, 直径2.0mm
用于电缆连接直径为2.0mm的根管治疗设备
E-Drive Apex Clip, pin diameter 2.0 mm
For endometric devices with a cable connection of Ø 2.0 mm

9940.000


E-Drive根尖夹, 直径1.5mm
用于电缆连接直径为Ø1.5mm的根管治疗设备
E-Drive Apex Clip, pin diameter 1.5 mm
For endometric devices with a cable connection of Ø 1.5 mm

9941.000

E-Drive喷雾器
用于清洁和润滑E-Drive根管治疗弯机
E-Drive Spray
For cleaning and lubricating the endodontic contra-angle E-Drive

LOCA 11.000

Locapex five根尖定位仪
在根管预备期间用于定位根尖的根管治疗测量设备
Locapex five
Endodontic measuring device for localizing the apex during the preparation of the root canal

Endo Rescue Kit – Removal of instrument fragments

The fracture of an instrument during an endodontic treatment not only causes the dentist enormous stress, it also poses an increased risk of post-endodontic complications to the patient. The removal of such fragments is often extremely difficult and almost impossible to plan.

The Endo Rescue Kit provides a simple and systematic solution, providing access to the opening of the root canal and allowing the removal of the fractured instrument. Once straight access to the fragment has been created by means of a conventional endodontic drill and two Gates burs, two specifically developed instruments greatly simplify a previously complicated procedure:

A centre drill (RKP) exposes the coronal part of the fragment. An extremely fine trepan bur is then placed onto the fragment which is seized by the bur and held in place by dentin residues. The fragment is then pulled out of the root in an anti-clockwise direction. Thanks to the small diameter of these two instruments, the fragment can be extracted whilst removing the absolute minimum of dentin.

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
9944


Clinical sequence

1. Initial situation. Presence of a fragment of a fractured instrument in a mesial root.
2. Recreation of the access cavity (H269GK.315.016)
3. The coronal curvature is straightened with axial movements to create direct access to the fragment (G180A.204.110).
4. Preparation of the access right down to the fragment (G180A.204.090)
5. Exposure of the fractured surface by drilling around the fragment (RKP.204.090)
6. The fragment is surrounded and seized. Anti-clockwise rotation. Removal of the fragment which is firmly held in the trepan bur by the residues of dentin (RKT.204.090)
7. The last third can now be correctly prepared.
8. Final situation after filling.

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
4601.000

Endo Rescue取根管断针套装
Endo Rescue Kit
For the removal of fractured instruments

	H269GK.315.016	1	
	G180A.204.110	1	
	G180.204.090	1	
	RKP.204.090	1	
	RKT.204.090	1	
	155.000.	1	

RKP


		2	2	2
尺寸 · Size	∅ $\frac{1}{100}$ mm	070	090	110

RA车针 (RA) RA · Right-angle (RA)

RKP.204. ... 070 090 110

∅_{max.} 20000 rpm
∅_{opt.} 300 rpm
Endo Rescue中心钻
Endo Rescue Centre Drill

RKT


		2	2	2
尺寸 · Size	∅ $\frac{1}{100}$ mm	070	090	110
D	∅ $\frac{1}{10}$ mm	7	9	11
D ₁	∅ $\frac{1}{10}$ mm	4	5	7

RA车针 (RA) RA · Right-angle (RA)

RKT.204. ... 070 090 110

∅_{max.} 20000 rpm
∅_{opt.} 300 rpm
Endo Rescue 环钻
逆时针旋转使用
Endo Rescue Trepan bur
To be used in anticlockwise rotation

Insert tray – antibacterial

The instruments of the AlphaKite and EasyShape systems and other endodontic instruments can be stored in a clearly arranged treatment tray. The insert tray is made of a special antimicrobial plastic.

- PP, 64 x 59 x 19 mm
- Antimicrobial effect thanks to the continuous release of silver ions
- High purity grade due to reduced bacterial contamination during idle times
- No formation of a biofilm on the surface of the material. No unpleasant smells.
- Reduced formation of breeding grounds for bacteria in hidden corners or depressions etc. Successfully used as a component in food processing machines, in the food and drink industry and in water technology.

Sterilization container made of high-performance plastic

- No seals, no maintenance required. More than 2,000 sterilization cycles
- Suitable for vacuum steam sterilization
- ePFTE filter. The long-lasting ePFTE filter remains in the sterilization container during reprocessing (manual or mechanical)
- Transparent material – the contents are visible from the outside
- Two or more containers can be laterally connected

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Universal Endo insert tray, for 28 endodontic instruments (handpiece and contra-angle), PP with silver ions, antibacterial (without instruments)

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Universal Endo kit, sterilisation container and insert tray (without instruments)

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Sterilisation filter 25 x 61 mm for sterilisation container A8, change after 12 months or after 150 sterilisation cycles, ePTFE, 2 pcs.

556.000.	1	A8灭菌容器 Sterilisation container A8
541.000.	1	通用型Endo插入盘 Universal Endo insert tray

9880

用于灭菌容器的日期插片，带有年份标记，在接近150次灭菌周期或至少一年后更换过滤器时调换或重新设定。日期插片没有以不同颜色显示
Date insert for sterilisation container, with indication of the year, to be exchanged or reset when the filter is changed after approx. 150 sterilisation cycles or at least once a year. The date insert comes in a different colour every year

Name	
Steril. Dat.	
Verf. Dat.	
Steril. Nr.	Dampf (braun)

9878

带点状标记的灭菌容器密封标签，可选配件，每次灭菌使用1个。灭菌期间点状标记改变颜色
Sealing label for sterilisation container with dot indicator, optional accessory, 1 label per sterilisation. The dot indicator changes colour during the sterilisation process


9879

灭菌容器的安全密封条，可选配件，每次灭菌使用1个
在盖子打开时，封条断开
Safety seal for sterilisation container, optional accessory, 1 seal per sterilisation
The seal breaks when the lid is opened

590 U.000

尺寸 · Dimensions mm 24,5 x 19 x 4,6

用于EasyShape插入盘的频率弹夹（例如：AlphaKite，通用的根管治疗插入盘）
用于记录根管锉和其他器械的使用次数
可以用标准的防水记号笔标记
Frequency clip for various insert trays (e.g. AlphaKite, universal endodontic trays)
To record how many times the endodontic files and other instruments have been used
Can be marked with standard waterproof permanent markers

590 ES.000

尺寸 · Dimensions mm 24,5 x 19 x 4,6

用于EasyShape插入盘的频率弹夹
用于记录EasyShape根管锉使用的次数
可以用标准防水记号笔标记
Frequency clip for EasyShape insert tray
To record how many times EasyShape files have been used
Can be marked with standard waterproof permanent markers

Wash box

The wash box is intended for mechanical cleaning and disinfection of instruments in the thermo disinfectant. Following an initial pre-cleaning step, the instruments to be cleaned are placed into the insert trays that are part of the chosen treatment system. The loaded insert tray can then be put into the wash box. It is then cleaned and disinfected in the thermo disinfectant.

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
Wash box

For mechanical cleaning and disinfection of instruments in the thermo disinfectant

new

595.000

60

带泡沫衬垫的根管治疗仪器中间支架 (5件)
用于治疗过程中中间卫生保存和清洁根管仪器 (无器械)
*Intermediate support for endodontic instruments with foam inserts (5 pcs)
For the hygienic intermediate storage and cleaning of root canal instruments during the treatment (without instruments)*

9866

尺寸 · Dimensions	mm	50 x 30 x 17
-----------------	----	--------------

白色泡沫衬垫, 25件替换装
Foam insert white, refill 25 pcs.

9890 L 4

尺寸 · Dimensions	mm	72 x 20 x 40
-----------------	----	--------------

作为通用仪器固定器的含有6个蓝色硅胶塞的通用型不锈钢车针架, 可以组合不同类型的柄, 用于最长为37mm的器械
Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 37 mm

9890 L 5

尺寸 · Dimensions	mm	72 x 20 x 50
-----------------	----	--------------

作为通用仪器固定器的含有6个蓝色硅胶塞的通用型不锈钢车针架, 可以组合不同类型的柄, 用于最长为47mm的器械
Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 47 mm

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Alpha Sequencer, instrument block made of stainless steel, teflon insert with 12 instrument holes (handle or right angle shank), 3 Sterimeters (yellow, red, blue) made of silicone for counting the sterilisation cycles already carried out, measurement scale on the inside of the lid

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Alpha Sequencer, instrument block made of stainless steel, teflon insert with 12 instrument holes (handle or right angle shank), 3 Sterimeters (green, black, white) made of silicone for counting the sterilisation cycles already carried out, measurement scale on the inside of the lid

9891


	8	8	8	8
尺寸 · Size	1	2	3	4
9891.000.	1	2	3	4

硅胶塞，含硅胶塞车针架的替换装
Silicone plug, refill for bur blocks with silicone plugs

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

Leakage test with dye

EasySeal (无染料渗透)
EasySeal (no dye penetration)

商用sealer X (无染料渗透)
Commercially available sealer X (no dye penetration)

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Natural tooth, prepared up to ISO 035, filled with different sealers and gutta-percha, immersed in dye penetrant for 48 hours, removed from the dye penetrant, cut in thin slices in intervals of 2 - 6 mm from the apex and examined for residues of the dye penetrant under a microscope.

EasySeal – Root filling material on the basis of epoxy resin

After a successful treatment of the root canal, the root filling is responsible for tightly sealing the canal to prevent reinfection.

A permanent apical seal of the root canal can be achieved with EasySeal, a root filling material based on epoxy resin. EasySeal is particularly suitable for orthograde filling techniques. It is self-sterile to prevent reinfection, radiopaque and dimensionally stable. EasySeal comes in a double-chamber syringe for easy, safe handling without need for previous mixing. One syringe contains 12 g which is enough for approx. 20 root fillings. The workable time of EasySeal is 15 minutes. Full setting is reached at 37°C body temperature within 24 hours.

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


RA车针 (RA) RA · Right-angle (RA)

GPXF21.204. ... 025 030

○_{max.} 4000 rpm
 无切缘的牙胶清除车针
 旋转产生的摩擦热导致牙胶塑性话，镍钛合金
 Gutta-percha remover without cutting edges
 Plastification of gutta-percha due to frictional heat caused by rotation, nickel-titanium alloy

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


Trepan bur for exposing fragments in the root canal

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


Composite System

DentinBuild Evo

is a dual curing composite for the cementation of root posts and for core build-ups. DentinBuild Evo is suitable for cementing and building up root posts made of glass fibre reinforced composite, ceramic and titanium.

DentinBond Evo

is a self-etching, dual curing adhesive for use with the DentinBuild composite.

DentinPost & Core Kit

To have all components ready at hand and to enable systematic work, the DentinPost & Core Kit contains everything required for post endodontic treatments in one case - ranging from the root post DentinPost Coated including all matching instruments for the preparation of the post site via the self-etching adhesive DentinBond Evo to the composite DentinBuild Evo for the cementation of root post and core build-ups, including useful core forms for core build-ups.

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
DentinPost & Core 套装 | 4613

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A large grid of small dots for writing, consisting of 20 columns and 30 rows of dots.

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MULTiflex® is a registered trademark of KaVo Dental GmbH, Germany
Brasseler DMI (Shanghai) Ltd.