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Brasseler®, Komet®, CeraBur®, CeraCut®, CeraDrill®,
CeraFusion®, CeraPost®, DC1®, DCTherm®, F360®,
F6 SkyTaper®, H4MC®, OptiPost®, PrepMarker®, PolyBur®,
TissueMaster®, TMC® i TissueMaster Concept® to
zastrzeżone znaki towarowe Gebr. Brasseler GmbH & Co. KG.

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*Brasseler®, Komet®, CeraBur®, CeraCut®, CeraDrill®,
CeraFusion®, CeraPost®, DC1®, DCTherm®, F360®,
F6 SkyTaper®, H4MC®, OptiPost®, PrepMarker®, PolyBur®,
TissueMaster®, TMC® und TissueMaster Concept® are
registered trademarks of Gebr. Brasseler GmbH & Co. KG.*

*Some of the products and designations mentioned in
the text are patented or copyrighted. The absence of a
special reference or the sign ® should not be interpreted as
the absence of copyright protection.*

*This publication is copyrighted. All rights, also with regard to
translation, reprint and reproduction (also in the form of
extracts) are reserved. No part of this publication may be
reproduced or reprocessed using electronic systems in any
form or by any means (photocopying, microfilm or other
methods) without the written permission of the editor.*

*Colours and products are subject to alterations.
Printing errors excepted.*

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General instructions for use and safety recommendations for the application of dental instruments, pins, posts and bur blocks

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Area of application

These general instructions for use and safety recommendations apply to all products and have to be generally observed. Non-observance of these instructions for use and safety recommendations increases the risk of injury and may impair the proper function of the instruments.

Separate instructions for use are enclosed in the packaging of products that require more detailed information (packaging is provided with the  symbol). These take precedence over the general instructions.

Storage

Protect packed instruments from UV rays and high temperatures. Store all instruments in a clean and dry environment. Do not store in the same room as solvents or chemicals.

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1. Proper use

- Those instruments that are supplied non-sterile have to be prepared prior to first use.
- Make sure that only technically and hygienically perfect and cleaned power systems (hand pieces and contra-angles) are used.
- Depending on their shank type, insert the instrument into the chuck as deeply as possible. Make sure that they are properly locked.
- The instrument must be rotating at the desired speed before contact is made with the tissue or material.
- Avoid jamming and using the instrument as a lever. Excessive contact pressure has to be avoided.
- For eye protection wear safety glasses. Use appropriate respiratory protection (mouth and nose). In the dental laboratory, use appropriate suction unit.
- Avoid contact with the instruments' working parts as this may increase the risk of injury.

Make sure to provide sufficient cooling by means of air/water spray. Additional external cooling is required when using instruments with extra-long shanks or oversized working parts. Damaged or corroded instruments have to be discarded.

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2. Recommended speeds

Please make sure to observe the recommendations for use and recommended speeds as indicated in the instructions of use and on the packaging of the products.

- $\bigcirc_{\max} 300\,000\text{ rpm}$ means: Suited for micro-motor contra-angles and turbines with stable ball bearings. Do not use in turbines with air bearing.
- $\bigcirc_{\max} 200\,000\text{ rpm}$ means: Suited for micro-motor hand pieces and contra-angles or lab hand pieces up to the speed indicated. Do not use in turbines.

Not observing the maximum permissible speed leads to an increased safety risk.

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7. Specific instructions for individual instrument types

- Avoid any contact with H_2O_2 (hydrogen peroxide). Make sure that the specified immersion times in the cleaning and disinfecting agents are not exceeded. The carbide working parts would be attacked (a typical indication is black staining of the instrument) reducing the instrument's service life.
- Tool steel instruments corrode and can therefore not be sterilized with a sterilization method using saturated steam (e. g. autoclave).
- To achieve an optimal surface roughness, subsequent finishing is necessary after using a diamond grinding instrument with coarse or very coarse grit.
- Use a disc guard for diamond discs when working intraorally.
- When using diamond separating strips and diamond strips please avoid contact with the gingiva as there is a risk of injury.
- Due to the ergonomically flat design of reciprocating files, please avoid jamming, bending or using the file as a lever during use. Please make sure that the files are properly fixed in the chuck of the reciprocating contra-angle to the required depth.
- To avoid overstressing of the instrument, root canal instruments made of nickel-titanium have to be used in a torque limited motor.
- Stainless steel root canal reamers (type “Gates”, “Müller”, “P”) are only intended for preparation of the coronal portion of the canal.
- For reprocessing anodized aluminium bur blocks, cleaning and disinfecting agents suitable for aluminium must be used. Other agents would destroy the anodized layer of the bur block. Aluminium blocks are generally not suited for preparation in the thermo disinfectant. Prior to sterilization, rinse bur block under running water and dry thoroughly (e.g. by air blasting). Frequent reprocessing may lead to colour deviations.
- Root canal posts made of fibre reinforced composite or PMMA cannot be sterilized. Therefore, they must be cleaned and disinfected with medical alcohol. Root posts are intended for single use only.
- Clean and sharpen DSB abrasives every now and then with the cleaning stone.
- Always use polishers and brushes with low pressure to minimize heat generation. Always polish in circular motion. To achieve brilliant high shine, use the polishers in the indicated sequence when using multiphase polishing systems.
- Polishers and Arkansas abrasives have to be prepared with suitable, alcohol-free agents (e.g. Komet DC1).


7.1. Instruments with internal cooling and sonic/ultrasonic instruments

- Manual reprocessing:
Rinse off surface contamination under running water. Remove all stubborn contamination with a nylon brush under water level, turning the instrument constantly. To assure patency of the internal cooling channel, penetrate the perforation with the cleaning wire. Rinse the cooling channel, e.g. with a syringe, until it is clean without any further contamination at all.
- Mechanical reprocessing:
Sonic and ultrasonic instruments are placed into the thermo disinfectant using the rinse adapter (see instructions for use).

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
8. Potential effects of using worn instruments

- *Discard any instruments with damaged or worn working parts as the use of damaged or worn instruments would have a negative effect on the work result.*
- *Damaged and deformed cutting blades will cause the instrument to vibrate and lead to poor preparation margins and rough surfaces.*
- *Void spots on the surface of diamond instruments are a sign of missing diamond particles and reduce the instrument's abrasive efficiency. An inferior diamond coating quality will result in excessive heat generation. Excessive contact pressure or temperatures may cause damage to the tissue.*
- *Please avoid excessive contact pressure during use as this may result in damage to the working part (Nicks on the blades, premature blunting and excessive heat generation).*
- *In abrasive instruments, excessive contact pressure can lead to stripping of the grit or clogging of the instruments and increased heat generation.*
- *To avoid undesirable heat generation during preparation, make sure to provide sufficient cooling by means of air/water spray (at least 50 ml/min).*
- *Additional external cooling is required when using instruments with a total length of more than 22 mm or a head diameter exceeding 2,5 mm.*
- *Not observing the maximum permissible speed will result in an increased risk of injury.*
- *There is an increased risk of infection in reusable instruments which have not been properly reprocessed.*
- *The reuse of disposable instruments (marked  on the packaging) is not permitted (e.g. polishers with lamellae and dental brushes). The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed (e.g. due to the risk of fracture with root canal instruments).*

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9. Safety and possible side effects

The above mentioned recommendations with respect to cooling, contact pressure, disinfection, cleaning and sterilization are to be strictly observed. The instruments should only be used for the intended application, as indicated by the symbols provided. Non-observance of these safety recommendations may lead to damage of the power system and/or injury, such as thermal necrosis, undesired preparation of tissue, damage to tissue or nerves, violation of the biological width, or infections. During preparation, some instruments may generate metallic abrasion which may lead to the presence of artifacts in MRI diagnostics.

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10. Liability

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

*Store products out of children's reach.
For dental use only.*

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¹ After the expiration date, the use of sterile packed instruments carries a risk of infection, and the proper function of the product can therefore not be guaranteed.

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² Do not use in case of damaged packaging as the sterility of the product is no longer guaranteed. Make sure to check that reusable instruments are in perfect condition and reprocess prior to use.

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³ The reuse of these products carries a risk of infection. A safe, risk-free use can therefore not be guaranteed.

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* With some instruments the grit size may deviate from the specified value, depending on their shape and size.

Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25 · 32657 Lemgo
Postfach 160 · 32631 Lemgo · Germany

Export:

Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-329
export@kometdental.de
www.kometdental.de