




F6 SkyTaper.
A new level of flexibility.



Brasseler®, Komet®, CeraBur®, CeraCut®,
CeraDrill®, GeraFusion®, CeraPost®, DC1®, DCTherm®, F360®,
F6 SkyTaper®, H4MC®, OptiPost®, PrepMarker®, PolyBur®,
TissueMaster®, TMC® and TissueMaster Concept® are
registered trademarks of the company Gebr. Brasseler GmbH &
Co. KG.

Some of the products and designations mentioned in the text
are patented or copyrighted. The absence of a special
reference or the sign ® should not be interpreted as the
absence of copyright protection.


This publication is copyrighted. All rights, also with regard to
translation, reprint and reproduction (also in the form of
extracts) are reserved. No part of this publication may be
reproduced or reprocessed using electronic systems in any
form or by any means (photocopying, microfilm or other
methods) without the written permission of the editor.

Colours and products are subject to alterations. Printing errors
excepted.

As at July 2018

Larger taper, limitless potential.

- Rotary preparation with one taper .06 file
- All-comprehensive range of files in 5 sizes to suit all root canal anatomies
- Flexibility and efficiency thanks to an intelligent file design
- Sterile packed single-use files



Flexibility without compromise.

Flexible system and application

- All-comprehensive system: NiTi files in five different sizes and 3 lengths to suit every canal anatomy
- Unrestricted choice of the obturation method thanks to the taper .06
- Rotary use in all torque-limited endodontic contra-angles and motors

Efficient and flexible file design



- One file for the majority of root canals
- Efficient double-S cross-section for thorough cleaning
- Excellent preservation of the course of the canal

[01] Generous preparation of the entrance of the root canal, e.g. with the opener OP10L19.


[02] Creation of a glide path, for example with a K-file in size 015.

The suitable F6 SkyTaper file is chosen according to the size of the previously used (manual) file.
Rule: Skip one ISO file size

[03] Preparation of the root canal with the F6 SkyTaper file.



File cross-section



Dr. Günther Stöckl

Specialist in endodontics,
Rottenburg an der Laaber, Germany.

“The F6 SkyTaper system is suitable for both general dentists and specialists. The excellent cutting performance of the instruments shortens the preparation time so that I can spend more time on the chemical disinfection.”

“I’m especially impressed with the 6% taper which prepares the canal systems particularly well for the subsequent hot or cold obturation. The F6 SkyTaper system has become an indispensable part of my preparation routine!””



F6 SkyTaper. Perfection in every detail.



Glide path up to #010

F06L21/L25/L31.204.020



Glide path up to #015

F06L21/L25/L31.204.025



Glide path up to #020

F06L21/L25/L31.204.030




+ Optional
for wide canals



F06L21/L25/L31.204.035



F06L21/L25/L31.204.040



Advantages of sterile packed single-use instruments:

- Reduced risk of fracture caused by cyclic fatigue
- The risk of cross-contamination is excluded
- No time spent on reprocessing (cleaning, disinfection and sterilization) and recording the frequency of use



9938 E-Drive
Endodontic contra-angle

Accessories for preparation:




● **OP10L19.204.030**
Opener



○ **PG03L21/L25.204.015**
PathGlider



● **PG03L21/L25.204.020**
PathGlider



595 Interim bur block

Accessories for obturation:



9994 F360 Fill Oven



9996 F360 Fill Obturators



9995 F360 Fill Verifier



PPF06 F6 SkyTaper paper points

GPF06 F6 SkyTaper gutta-percha points



9978 EasySeal

Root filling material based on epoxy resin

Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25 · 32657 Lemgo
Postfach 160 · 32631 Lemgo · Germany

Verkauf Deutschland:
Telefon +49 (0) 5261 701-700
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

Export:
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-329
export@kometdental.de
www.kometdental.de
