# EnGO

## Endo Motor

EN\_Instruction Manual

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Brasseler<sup>®</sup>, Komet<sup>®</sup>, Art2<sup>®</sup>, CeraBur<sup>®</sup>, CeraCut<sup>®</sup>, CeraDrill<sup>®</sup>, CeraFusion<sup>®</sup>, CeraPost<sup>®</sup>, DC1<sup>®</sup>, DCTherm<sup>®</sup>, Derminator<sup>®</sup>, FastFile<sup>®</sup>, F360<sup>®</sup>, F6 SkyTaper<sup>®</sup>, H4MC<sup>®</sup>, OccluShaper<sup>®</sup>, OptiPost<sup>®</sup>, PolyBur<sup>®</sup>, PrepMarker<sup>®</sup>, Procodile<sup>®</sup>, R6 ReziFlow<sup>®</sup>, TissueMaster<sup>®</sup>, TMC<sup>®</sup>, TissueMaster Concept<sup>®</sup> and Visio-Soft<sup>®</sup> are registered trademarks of the company Gebr. Brasseler GmbH & Co. KG.

ReFlex<sup>®</sup> und EndoPilot<sup>®</sup> are registered trademarks of the company Schlumbohm GmbH & Co. KG.

Some of the products and designations mentioned on our web sites are protected by patent law, trademark law and copyright. The absence of an express statement or of the symbols<sup>®</sup> should not be interpreted as absence of protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or reprocessed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products are subject to alterations.

Printing errors excepted.

As at January 2023

# Introduction Manual

1 Product introduction .....	1
2 Installation .....	2
3 Function and operation of product .....	5
4 Operation instruction .....	7
5 Troubleshooting .....	13
6 Cleaning, Disinfection and Sterilization .....	13
7 Storage, maintenance and transport .....	17
8 Environmental protection .....	17
9 After sales service .....	17
10 European authorized representative .....	17
11 Icons .....	17
12 Statement .....	18
13 EMC-Declaration of conformity .....	18

# 1 Product introduction

## 1.1 Product description

The endo motor (EnGO motor) is mainly used in endodontic treatments. It is a cordless endo motor with an incorporated root canal measuring function. It can be used as an endo motor for the preparation and enlargement of root canals or as a device for measuring the canal length. The EnGO motor is capable of monitoring the position of the file tip inside the canal while enlarging the canal.

Features:

- a) Efficient brushless motor, low noise, long service life.
- b) Cordless portable endo motor with incorporated length determination.
- c) 360 degree rotation of the contra angle.
- d) Make use of real-time feedback technology and dynamic torque control, to effectively prevent inadvertent file separation.

## 1.2 Model and specification

Model: Ai-Motor(EnGo-Motor)

Please refer to packing list for device configurations.

## 1.3 Components

The device is composed of charging base, motor handpiece, contra angle, measuring wire, lip hook, file clip, power adapter, protective silicon cover, etc.

## 1.4 Indications for Use

The EnGO motor is a cordless motorized handpiece for endodontic treatments with an incorporated root canal measuring feature. It can be used for the preparation and enlargement of root canals. The EnGO motor is capable of monitoring the position of the file tip inside the canal while enlarging the canal.

## 1.5 Scope of application

The device may only be operated by qualified dentists in hospitals and clinics.

## 1.6 Caution

Federal law restricts this device to sale by or on the order of a dentist.

## 1.7 Contraindications

- a) Not to be used by doctors with cardiac pacemakers.
- b) The device must not be used on patients with hemophilia.
- c) Use with caution in patients with heart disease, pregnant women and young children.

## 1.8 Warnings

- 1.8.1 Please read this Instruction Manual carefully before first operation.
- 1.8.2 This device should be operated by professional and qualified dentists in a qualified hospital or clinic.
- 1.8.3 Do not directly or indirectly place this device near heat source. Operate and store this device in a safe environment.
- 1.8.4 This device requires special precautions regarding electromagnetic compatibility (EMC) and must be in strict accordance with the EMC information for installation and use. In particular, do not use this equipment in the vicinity of fluorescent lamps, radio transmitting devices, remote control devices, handheld and mobile high-frequency communication devices.
- 1.8.5 Please use the original contra angle. The use of other contra angles is not authorized and may have adverse consequences.
- 1.8.6 Please do not make any changes to the device. Any changes may violate safety regulations and cause harm to the patient. No liability will be accepted in cases where changes were made to the device.
- 1.8.7 Please use original power adapter. The use of other power adapters will result in damage to the lithium battery and control circuit.
- 1.8.8 The motor handpiece cannot be autoclaved. Use disinfectant of neutral pH value or ethyl alcohol to wipe its surface.
- 1.8.9 Do not press the push cover of the contra angle until this has stopped rotating. Failure to observe this will result in the contra angle breaking.
- 1.8.10 Do not remove the contra angle until the motor handpiece has stopped rotating. Failure to observe this will result in the contra angle and the gear inside the motor handpiece breaking.
- 1.8.11 Please ensure that the file is properly positioned and locked before starting the motor handpiece.
- 1.8.12 Please set torque and speed as per the recommended specifications of file manufacturer.
- 1.8.13 Incorrect replacing of lithium batteries can lead to unacceptable risks, so use the original lithium battery and replace the lithium battery according to the correct steps in the instructions.
- 1.8.14 Please remove the battery if the motor handpiece is not likely to be used for some time.
- 1.8.15 Wireless charging will generate heat, and the surface temperature of the charging base and the motor handpiece will rise. It is recommended that the time of contact between the motor handpiece and any charging base during wireless charging should not exceed 10 seconds (only for EnGOMotor).

## 1.9 Device safety classification

- 1.9.1 Type of operation mode: Continuously operating device
- 1.9.2 Type of protection against electric shock: Class II equipment with internal powered electronic medical supply
- 1.9.3 Degree of protection against electric shock: B type active part
- 1.9.4 Degree of protection against harmful ingress of water: Ordinary equipment (IPX0)

1.9.5 Degree of protection during use in the presence of flammable mixtures of anaesthetics with air, oxygen or nitrous oxide: The device must not be used in the presence of flammable mixtures of anaesthetics with air, oxygen or nitrous oxide.

1.9.6 Active parts: contra angle, lip hook, file clip, touch probe.

1.9.7 The contact duration of active parts: 1 to 10 minutes.

1.9.8 Max. temperature of the surface of active parts may reach 46.6°C.

## 1.10 Primary technical specifications

### 1.10.1 Battery

Lithium battery in motor handpiece: 3.7V /2000mAh

### 1.10.2 Power adapter

(Model: ADS-6AM-06N 05050/UE08WCP- 050100SPA)

Input: ~100V-240V 50Hz/60Hz 0.4A Max

Max Output: 5.0V1A

1.10.3 Torque range: 0.4Ncm-5.0Ncm(4mNm~50mNm)

1.10.4 Speed range: 100 rpm ~ 2500 rpm

1.10.5 Wireless charging

Frequency range: 112-205KHz


Maximum RF output power of the product: 9.46dBuA/m@3m

## 1.11 Environment parameters

1.11.1 Ambient temperature: +5°C ~ +40°C


1.11.2 Relative humidity: 30% ~ 75%

1.11.3 Atmospheric pressure: 70kPa ~ 106kPa


## 2 Installation


### 2.1 Basic accessories of product


Motor handpiece


Contra angle


EnGO Motor Charging base


Nozzle

Measuring wire


Protective silicon cover


Power adapter

Touch probe

Disposable insulation sleeves

File clip

Lip hook

### 2.2 Display Screens

2.2.1 Display Screens for 6 Operation Modes and Standby

2.2.1.1 EAL Mode


This mode is for canal measurement.  
The motor handpiece does not run in this mode.

2.2.1.2 CW Mode

The motor handpiece rotates forward 360°, clockwise direction.


2.2.1.3 CCW Mode

The motor handpiece rotates in a counterclockwise direction only.  
This mode is used to inject calcium hydroxide or other medicaments.  
When this mode is used, a double-beep sounds continuously.


2.2.1.4 RECI Mode Safety Glide Path Mode

F: Forward angle, R: Reverse angle  
The rotation angle is adjustable, but the forward angle must be equal to the reverse angle.

2.2.1.5 ATR Mode


ATR: Adaptive Torque Reverse function.

Normal continuous forward rotation, when the load of the file is greater than the set torque limit, the file will start to rotate alternately at the set angle.


2.2.2 Torque Display

This appears when the motor is running.  
Meter shows the torque load on the file.

2.2.3 Canal Measurement Display

This appears when a file is inside the canal and the lip hook is contacting the patient mouth. The bars in the meter show the location of the file tip.  
In the EAL Mode, if the length is less than 1.0, the display will be enlarged.

The meter numbers 1.0, 2.0, 3.0 and digital numbers 00 -16 do not represent the actual length from the apical foramen. They simply indicate the

file progression towards the apex. The digital numbers -1 and -2 indicate that the file has passed the apex foramen.

The digital number "00" indicates that the file has reached the apex foramen.

Subtract 0.5 -1 mm from the measured file length to calculate the working length. These numbers are used to estimate the canal's working length.

## 2.3 Instructions for contra angle

2.3.1 The contra angle uses precision gear transmission, and the transmission ratio is 6:1.

2.3.2 Before the first use and after treatments, please clean and disinfect contra angle with disinfectant of neutral PH value. After disinfection, lubricate it with specific cleaning oil. Finally, sterilize it under high temperature and (0.20 MPa ~ 0.23 MPa).

2.3.3 The contra angle can only be used in combination with this device. Otherwise the contra angle will be damaged.


## 2.4 Installation and removal of contra angle

### 2.4.1 Installation

Align any locating pin of the contraangle with the positioning slot on the motor handpiece and push the contra-angle horizontally. The three locating pins on the contra-angle are inserted into the three positioning holes on the motor handpiece.

A "click" sound indicates that the installation is in place. The contraangle can be freely rotated by 360°.

The contra-angle is free to rotate, adapting to the root canal in different positions, and it is advisable to watch the screen when operating.

### 2.4.2 Removal

Pull out the contra angle horizontally when the motor handpiece is not in operation.

Warnings:


- Before plugging in or pulling out contra angle, please first stop the motor handpiece.
- After installation, please check and confirm that the contra angle is firmly locked in place.

## 2.5 Installation and removal of file

### 2.5.1 Installation of file

Before starting the device, plug the file into the hole of the contra angle head.

Hold down the push button on the contra angle and insert the file. Turn the file back and forth until it is lined up with interior locking groove and slips into place. Release the button to allow the file to engage in the contra angle.

Warnings:

After plugging the file into the contra angle, guide your hand across the push cover to assure that the file cannot be taken out. Be careful when inserting files to avoid injury to fingers.


Inserting and removing files without holding down the push button may damage the chuck of the contra angle.

Please use files with shanks that meet the ISO standard.

(ISO standard:  $\varnothing 2.334 - 2.350$  mm)

### 2.5.2 Removal of file

Press down the push cover, and then directly pull out the file.

Warnings:

Before plugging in and pulling out the file, the motor handpiece must be stopped. Be careful when removing files to avoid injury to fingers. Removing files without holding down the push button will damage the chuck of contra angle.

## 2.6 Connection of the canal measurement equipment


This is not required if the canal measurement function will not be used.

Connect the measuring wire to the motor handpiece. Line up the measuring wire plug with the notch on the back of the motor and push it all the way in. Connect the file clip plug with the socket (black) on the measuring wire. Connect the lip hook to the socket (white) on the measuring wire.

Warnings:

Connect the lip hook to the socket (white) on the measuring wire.

Otherwise, the functions of root canal preparation and root canal length measurement cannot be used together.

## 2.7 Installation and removal of disposable insulation sleeves

### 2.7.1 Installation

Before each use of the handpiece and after the handpiece is cleaned and disinfected, put on a disposable insulation sleeve. Take the insulation sleeve out of the insulation sleeve box, then insert the insulation sleeve into the motor handpiece from the thin end of the handpiece, and adjust the fit of the insulation sleeve until there are no more obvious wrinkles.

After installing the disposable insulation sleeve, wrap the barrier film around the handpiece surface. After that, clean and disinfect the surface of the handpiece. Refer to Chapter 6.3 for cleaning and disinfection procedures.


### 2.7.2 Removing

After each use, remove the barrier film and slowly pull off the insulation sleeve from the thin end of the handpiece.

Warning:

Insulation sleeves are not reusable.

## 3 Function and operation of product

### 3.1 Button definition and settings

#### a) Turn power on

Press main button to turn on motor handpiece.

#### b) Turn power off

Hold down the setting button "P", then press main button to turn off motor handpiece.

#### c) Customized program change

Press adjusting button "+" / "-" during standby mode.

#### d) Parameter setting

Press setting button "P" to select the required target parameters, press adjusting button "+" / "-" to change, then press main button or wait 5 seconds to confirm.

#### e) Preset program selection

Press setting button "P" for a longer period of time to enter the preset program from standby mode. Press the adjusting button "+" / "-" to select the required file system, then press the setting button "P" to confirm.

Next, press the adjusting button "+" / "-" to select a file or the number of files, then press the main button to confirm.

#### f) Handpiece functions setting


With the motor handpiece turned off, hold down the setting button “P” and press main button to enter handpiece functions setting, press setting button “P” until the required target setting has been reached, press adjusting button “+”/“-” to adjust, then press main button to confirm.

### 3.2 Screen display

Standby interface

- a) Customized program sequence number 0 - 9, totally 10 programs in total.
- b) Battery consumption
- c) Set speed
- d) Set torque
- e) Operation mode


Working interface

- a) Set speed
- b) Set torque
- c) Real time torque
- d) Torque display scale

Canal measurement mode interface

- a) Apical reference point flash bar
- b) EAL: Electronic apex locator


Canal measurement mode interface

- a) Canal length indicator bar
- b) Indication number

Digital numbers 00-16 do not represent the actual length from the apical foramen.

They simply indicate the file progression towards the apex. Number “00” indicates that the file has reached the apical foramen.

- c) Apical foramen.

Interface for setting the apical reference point

- a) Apical reference point flash bar
- b) Apical foramen
- c) Digital “02” meter reading, very near physiological apical foramen.

### 3.3 Terms and definitions

CW	Clockwise rotation, forward rotation Be applied to rotary file
CCW	Counter clockwise rotation, reverse rotation Be applied to special file, inject calcium hydroxide and other solutions
RECIP	Safety Glide Path Mode
ATR	Adaptive torque reverse Up to setting torque, the motor will move with ATR mode; when torque reduce to normal value, the motor will clockwise rotate.
Forward Angle	Angle of clockwise rotation of the file.
Reverse Angle	Angle of counter clockwise rotation of the file.
EAL	Electronic apex locator In the mode, the device will work like a stand-alone apex locator
AP	Apical foramen.
Apical Action	The file action when file tip reaches the flash bar point.
Flash Bar Position	Shows the point inside the canal where specified apical action is triggered.
Auto Start	The file rotation starts automatically when the file is inserted in the canal.
Auto Stop	The file rotation stops automatically when the file is taken out of the canal.
Apical Slow Down	The file slows down automatically as it approaches the apex. Activating in CW and CCW operation mode.
Operation Mode	5 operation modes for canal shaping and measurement. Such as CW, CCW, RECIP, ATR and EAL.

Speed	File rotation speed.
Torque (Torque Limit/Trigger Torque)	For CW and CCW modes, the torque value (Torque Limit) that triggers reverse rotation. For ATR mode, the torque value (Trigger Torque) that triggers ATR action.

## 4 Operation instruction

### 4.1 Power on and power off

#### 4.1.1 Starting and stopping of motor handpiece


a) While the handpiece is set to "Power off", press the main button and the handpiece will switch to standby mode. The display on the screen is as follows:

Standby mode

b) While in standby mode, press the main button and the handpiece will switch to operational mode. The display on the screen is as follows:

Operational mode

c) Press the main button again, and then the motor handpiece switches back to stand-by mode.

d) Hold down the setting button "P", then press main button to turn off motor handpiece.

In standby mode, the motor handpiece would automatically shut down after 3 minutes. The motor handpiece will also automatically shut down while it is put into the charging base.

### 4.2 Selecting customized program sequence number

The motor handpiece has 10 memory programs (M0-M9) and 5 preset programs, press adjusting button "+" / "-" to change customized program sequence number during standby state. M0-M9 are memory programs for canal shaping and measurement, every memory program has its own parameters such as operation mode, speed and torque, all these parameters can be changed.

### 4.3 Parameter setting

	<p>Before starting of motor handpiece, please check the operation mode is correct.</p> <p>All the parameters must be set according to files, make sure all the parameters are excepted before starting of motor handpiece, otherwise has risk of file separate.</p>
	<p>It has 5 operation modes for canal shaping and measurement: CW, CCW, RECIP, ATR and EAL (See chapter 3.3 Terms and definition for the explanations of these modes.) Press setting button "P" once during standby state, press adjusting button "+" / "-" to select correct operation mode. CCW mode is used to inject calcium hydroxide or other medicaments.</p> <p>When this mode is being used, a doublebeep sounds continuously, used to indicate that the device is rotating in a counter clockwise direction.</p>
<p>Repeatedly press the setting button "P" to check if all parameters of this operational mode are as required. If this is not the case, press the adjusting button "+" / "-" to select the correct parameters.</p>	
	<p>The speed setting can be adjusted from 100 rpm to 2500 rpm. Press adjusting button "+" / "-" to increase or decrease speed. Keep button pressed over a longer period of time to increase or decrease the speed fast.</p> <p>In ATR mode, speed of 100 ~ 500 rpm are available.</p> <p>In RECIP mode, speed of 100 ~ 500 rpm are available.</p>
	<p>The torque setting can be adjusted from 0.4 Ncm to 5.0 Ncm.</p> <p>Press adjusting button "+" / "-" to increase or decrease torque. Keep button pressed over a longer period of time to increase or decrease the torque fast. In ATR mode, the Trigger Torque of 0.4 Ncm ~ 4.0 Ncm are available.</p> <p>In RECIP mode, the torque of 2.0 Ncm ~ 5.0 Ncm are available.</p>

<p>Apical Action</p> <p><b>OFF</b></p>	<p>Actions that happen automatically when the file tip reaches the predetermined apical reference point inside the canal (flash bar position). Benefit from integration of length determination: when the file reaches the reference point, the motor will respond according to setting, i.e. go into reverse mode, stop or off. Press adjusting button “+”/“-” to change. OFF: To disable the apical action function. The file keeps rotating as before if when the reference point has been reached. Stop: The file automatically stops rotating when the reference point is reached. It will start rotating again if it is moved upward a little. Reverse: The file automatically switches to reverse mode when the reference point is reached or has been passed. It will change back to the previous rotational direction if it is moved upward a little.</p>
<p>Auto Start</p> <p><b>OFF</b></p>	<p>Rotation starts automatically when the file is inserted into the canal and the canal length indicator bar lights up more than 2 bars. Press adjusting button “+”/“-” to change. OFF: Motor does not start when file is inserted into the canal. The main button is used to start and stop the motor handpiece. ON: Motor starts automatically.</p>
<p>Auto Stop</p> <p><b>OFF</b></p>	<p>Rotation stops automatically when the file is taken out of the canal and the canal length indicator bar lights up less than 2 bars before the file is taken out. Press adjusting button “+”/“-” to change. OFF: Motor does not stop when file is taken out the canal. The main button is used to start and stop the motor handpiece. ON: Motor stops automatically.</p>
<p>Flash Bar Position</p> 	<p>This is the reference point where various apical actions are triggered. Press adjusting button “+”/“-” to select reference point by changing the flash bar. The meter's 0.5 reading indicates that the file tip is located very near the physiological apical foramen. The reference point (flash bar) can be set from 2 to AP (Apex) on the meter.</p>
<p>Apical Slow Down</p> <p><b>OFF</b></p>	<p>Rotation automatically slows down as the file tip approaches the reference point. Press adjusting button “+”/“-” to change. OFF: Disable Apical Slow Down function. ON: Rotation automatically slows down as the file tip approaches the reference point.</p>
<p>Forward Angle</p> <p><b>30°</b></p> <p>Reverse Angle</p> <p><b>30°</b></p> <p>M1 F:30°</p> <p>RECIP R:30°</p>	<p>Forward Angle In the RECIP mode, the forward angles of 20° ~ 400° are available. In the ATR mode, forward angles of 60° ~ 400° are available. Reverse Angle In the RECIP mode, reverse angles of 20° ~ 400° are available. In the ATR mode, reverse angle cannot be greater than the forward angle.</p>

#### 4.4 Preset program selection

<p>FQ06 300rpm  <input type="checkbox"/> CW 2.4Ncm</p>	<p>For convenience, we preset some common file systems.          Press adjusting button “+”/“-” to switch to preset program (M0-M9, preset program 1 - 5 ), the interface will be displayed as shown on the left.</p>
<p>FQ08 Opener          FQ03 Pathglider          FQ04 020-035          FQ04 045-055</p>	<p>Keep the setting button “P” pressed for a longer period of time “P” to enter preset program during standby mode, the interface will be displayed as shown on the left. Press adjusting button “+”/“-” to select file system.          The parameters of “FQ04” can also be altered to change them from the default setting.</p>
<p>FQ04 020-035          020-035 CW          300rpm          1.8Ncm</p>	<p>If you want to change back to default setting, keep the setting button “P” pressed down for a longer period of time to enter the preset program during standby state, select “FQ06” and press “Main” button to confirm, the default setting will be reloaded.          The default settings of the preset program can also be restored by turning the handpiece off and then back on.          Changing the preset program default setting is not recommended.</p>


#### 4.5 Handpiece functions setting

With the motor handpiece turned off, hold down the setting button “P” and press main button to enter handpiece functions setting, press setting button “P” till target setting, press adjusting button “+”/“-” to adjust, then press main button to confirm.

<p>Software Version  <b>V1.0.0</b></p>	<p>With the motor handpiece turned off, hold down the setting button “P” and press main button to enter handpiece functions setting, the software version number will appear on the display screen.</p>
<p>Auto Power OFF  <b>5 min</b></p>	<p>After 3 seconds of displaying the version number on the screen, the “Auto Power OFF” can be change, press adjusting button “+”/“-” to adjust, then press to “Main” button to confirm.          The device is equipped with an automatic switch-off function without need to press any buttons. This can be set from 3 to 30 minutes in 1-minute increments.</p>
<p>Auto Standby Scr  <b>30 sec</b></p>	<p>Press setting button “P” again, the “Auto Standby Scr” can be change, press adjusting button “+”/“-” to adjust, then press to “Main” button to confirm.          The handpiece motor automatically returns to standby mode without need to press any buttons.          This can be set from 3 to 30 seconds in 1-second increments.</p>
<p>Dominant Hand  <b>Right</b></p>	<p>Press setting button “P” again, the “Dominant Hand” can be changed, press adjusting button “+”/“-” to adjust, then press to “Main” button to confirm.          The right hand and the left hand can be set.</p>
<p>Calibration  <b>OFF</b></p>	<p>Press setting button “P” again, the “Calibration” can be change, press adjusting button “+”/“-” to select “ON”, then press to “Main” button to calibration.          Before calibrating, making sure the original contra angle is installed, and do not insert any files. The torque will not be correct if the calibration is carried out without the original contra angle and if the contra angle chuck is subject to any loads, and has risk of file separate. After replacement of contra angle, the contra angle shall be calibrated before use.</p>
<p>Beeper Volume  <b>Vol.3</b></p>	<p>Press setting button “P” again, the “Beeper Volume” can be changed, press adjusting button “+”/“-” to adjust, then press to “Main” button to confirm.          The “Beeper Volume” can be set from 0 - 3.          Vol.0: Mute.</p>
<p>Restore Defaults  <b>OFF</b></p>	<p>Press setting button “P” again, the “Restore Defaults” can be changed, press adjusting button “+”/“-” to select “ON”, then press to “Main” button to restore defaults.</p>

#### 4.6 Protective function of automatic reverse

During operation, if the load value exceeds the preset torque value, the file rotational mode will automatically change to reverse Mode. The file will once again return to normal rotational mode when the load is below the preset torque value again.

**Cautions:**

1. Protective function of automatic reverse is **ONLY** suitable for CW mode.
2. This function is forbidden under CCW mode, ATR mode.
3. When the motor handpiece battery indicator indicates a low battery capacity, the low battery capacity is insufficient to support the motor handpiece enough to reach the limit torque value, that is, the auto-reverse function will not work properly. Please charge it in time.
4. If the motor handpiece subjected to continuous loads, the machine may stop automatically as a result of overheat protection. If it happens, turn off the motor handpiece for a while until the temperature drops.

**4.7 Motor operation**

Please set operation mode, torque and speed as per the recommended specifications of file manufacturer.

	<p><b>Motor alone mode</b>                  When the device is used as a motor only, the torque bar will show on the screen.                  (for more information about torque bar, please see chapter 3.2 Screen display)</p>
	<p><b>Motor combined with canal measurement function mode</b>                  When using the motor in combination with the canal measurement function, the measuring wire must be connected with motor handpiece by USB socket, and the socket must be connected with the patient's lip by lip hook. The black socket is kept idle. The canal length indicator bar will show on the screen (for more information about canal length indicator bar, please see chapter 3.2 Screen display) Set parameters of automatic functions as needed, such as Apical Action, Auto Start, etc (more information about automatic functions, please see chapter 4.3 Parameter setting).</p>
	<p><b>Connection testing</b>                  It is strongly recommended to test the connection before every use. Touch the lip hook with the file in the contra angle and check that all the bars on the meter on the screen light up, and the motor should be reversed continuously, otherwise, the measuring wire or contra angle should be replaced.</p>

**4.8 Canal measurement procedure**

When the device is used as an apex locator only, we suggest that you put the motor handpiece on the charging base to get a better visual angle.  
 Press setting button "P" once during standby state, press adjusting button "+" / "-" to select EAL operation mode, then press main button to confirm. (See chapter 3.3 Terms and definition for information on the operation modes.) The measuring wire must be connected with motor handpiece by USB socket, socket is connected with the patient's lip by lip hook, and black socket is connected with file clip.  
 The canal length indicator bar will show on the screen (for more information about canal length indicator bar, please see chapter 3.2 Screen display).

The file clip must hold the file correctly.  
 Push the button on the file clip with your thumb in the direction shown by the arrow. Clip the holder onto the metal upper part of the file and then release the button.

OK NG

Connection testing  
 It is strongly recommended to test the connection before every use. Clip the holder onto lip hook and check that all the bars on the meter on the screen light up, otherwise, the measuring wire or file clip should be replaced.

**Root canals not suitable for canal measurement**  
 An accurate measurement of the root canal cannot be achieved under the below described conditions:

Root canal with a large apical foramen  
 Root canals that have an exceptionally large apical foramen due to a lesion or incomplete development cannot be accurately measured.  
 The result obtained by the measurement may be shorter.

Root canal with blood leaking from the opening  
 If blood oozes out of the opening of the root canal and comes in contact with the gums, this will result in electrical leakage and an accurate measurement cannot be obtained. Wait for bleeding to stop completely. Clean the inside and opening of the canal thoroughly to get rid of all blood, and then make a measurement.  
 Root canal with a chemical solution leaking from the opening.  
 An accurate measurement cannot be obtained if some chemical solution is leaking from the canal opening. In this case, clean the canal and its opening.  
 It is important to get rid of any solution leaking the opening.

Broken crown  
 If the crown is broken and a section of the gingival tissue intrudes into the cavity surrounding the canal opening, contact between the gingival tissue and the file will result in electrical leakage and an accurate measurement cannot be obtained. In this case, build up the tooth with a suitable material to insulate the gingival tissue.

gypsum

	<p><b>Fractured tooth</b> A fractured tooth will cause electrical leakage and an accurate measurement cannot be obtained. A branch canal will also cause electrical leakage.</p>
	<p><b>Re-treatment of a root filled with gutta-percha</b> The gutta-percha must be completely removed to eliminate its insulating effect. After removing the gutta-percha, pass a small file all the way through the apical foramen and then put a little saline solution in the canal, but do not let it overflow the canal opening.</p>
	<p><b>Crown or metal prosthesis touching gingival tissue</b> Accurate measurement cannot be obtained if the file touches a metal prosthesis that is touching gingival tissue. In this case, widen the opening at the top of the crown so that the file will not touch the metal prosthesis before taking a measurement.</p>
	<p><b>Extremely dry canal</b> If the canal is extremely dry, the meter may not move until it is quite close to the apex. In this case, try moistening the canal with saline solution.</p>
<p><b>Differing measuring results between apex locator reading and radiography.</b> Sometimes the reading of apex locator and the X-ray image will not correspond. This does not mean that the apex locator is not working properly or that the X-ray exposure is a failure. An X-ray image might not show the apex correctly depending on the angle of the X-ray beam, and the location of the apex might seem to be different than it really is.</p>	
	<p>The actual apex of the canal is not the same as the anatomical apex. There are frequently cases where the apical foramen is located up towards the crown. In these cases, an X-ray might indicate that the file has not reached the apex even though it has actually reached the apical foramen.</p>

#### 4.9 Battery Charging


The motor handpiece has built-in rechargeable lithium battery. When charging the battery, leave approximately 10cm around the charging base for easy access to inlet and the power cord. Insert the power adapter plug into the charging base power socket and confirm that they are correctly connected. Then insert the motor handpiece into the charging base (the motor handpiece needs to be correctly aligned with the charging base in the same direction for charging). When the blue indicator on the charging base flashes, it is charging. When the motor handpiece is fully charged, the blue indicator on the charging base remains permanently lit without flashing (Only for EnGO-Motor). After charging, please unplug the power adapter.

#### 4.10 Replacing Battery

Replace the battery if it seems to be running out of power sooner than it should. Please use the original lithium battery.

- a) Turn the motor handpiece power off.
- b) Use tweezers etc. to open the rubber cover and then remove the screw.
- c) Remove the battery cover.
- d) Remove the old battery and disconnect the connector.
- e) Connect the new battery and put it in the motor handpiece.
- f) Replace the cover and its screw.

It is recommended to contact local distributors or manufacturers to replace the battery.

#### 4.11 Oiling of contra angle

Only the original oil injection nozzle can be used for oiling of contra angle.

The contra angle needs to be lubricated after cleaning and disinfection, but before sterilization.

1. Firstly, screw the injecting nozzle onto the jet of oil bottle. (Around 1 to 3 turns)

2. Next, plug the nozzle into the end part of contra angle, and then grease the contra angle for 2-3s till the oil leaks out of the contra angle head part.

3. Place the end part of the contra angle in a vertical position for at least 30 minutes to allow any excess oil to escape under gravity.

Warnings: ⚠

Motor handpiece cannot be filled with oil.

Cautions:

a) To prevent the contra angle from detaching itself due to pressure, use your hand to hold the contra angle in place during lubrication.

b) Do not use a swirling nozzle. Swing nozzle can only be used for injection of gas, not for oiling.

## 5 Troubleshooting

### 5.1 Troubleshooting

Failure	Possible cause	Solutions
The motor handpiece does not rotate.	The device is in EAL mode, EAL mode is only for canal measurement.	Changing to CW, CCW, RECIP or ATR mode.
There is a continuous beeping sound after starting the motor handpiece.	The continuous beeping sound is indicating that the motor handpiece is in CCW mode.	Stop the motor handpiece and change the operating mode to CW Mode.
Contra angle calibration failure	Calibration failure caused by strong resistance of contra angle	Clean the contra angle, and recalibrate after oil injection.
The time of operation becomes shorter after charging.	Decreasing battery capacity.	Please contact local distributor or manufacturer.
No sound	Beeper Volume set to 0. Vol.0: Mute.	Set Beeper Volume to 1,2,3.
The continuously rotating file is stuck at the root canal.	Incorrect specification setting. The file is subjected to excessive torque loads.	Choose CCW Mode, start the motor handpiece, and take the file out.

## 6 Cleaning, Disinfection and Sterilization

### 6.1 Foreword

For hygiene and sanitary safety purposes, the contra-angle, the lip hook, the file clip, the protective silicon cover and the touch probe must be cleaned, disinfected and sterilized before each usage to prevent any contamination.

This concerns the first use, as well as all subsequent uses.

### 6.2 General recommendations

6.2.1 Use only a disinfecting solution which is approved for its efficacy (VAH/ DGHM-listing, CE marking, FDA and Health Canada approval) and in accordance with the DFU of the disinfecting solution manufacturer.

6.2.2 Do not place the contra-angle in a disinfectant solution or in an ultrasonic bath. Do not use chloride detergent materials.

6.2.3 Do not use bleach or chloride disinfectant materials.

6.2.4 For your own safety, please wear personal protective equipment (gloves, glasses, mask).

6.2.5 The user is responsible for the sterility of the product for the first cycle and each further usage. The user is equally responsible for ensuring that no damaged or contaminated instruments are used after sterilization, if applicable.

6.2.6 The water quality has to comply with the local regulations especially for the last rinsing step or with a washer-disinfector.

6.2.7 To sterilize the endodontic files, refer to the manufacturer's instructions for use.

6.2.8 The contra-angle needs to be lubricated after cleaning and disinfection, but before sterilization.

### 6.3 Cleaning and disinfection steps for the motor handpiece, the AC adapter and the base

Before and after each use, all the objects that were in contact with infectious agents should be cleaned using wipes impregnated with a disinfecting and detergent solution (a bactericidal, fungicidal and aldehyde free solution) approved by VAH/DGHM listing, CE marking, FDA and Health Canada solution manufacturer.

Warning:

Do not sterilize the motor handpiece, the AC adapter and the base.

#### 6.3.1 Pre-Op processing

Before each use, the handpiece, charger, and base must be cleaned and disinfected. The specific steps are as follows:

Warning:

The handpiece, charger, and base cannot be cleaned and disinfected with automatic equipment. Manual cleaning and disinfection is required.

##### 6.3.1.1 Manual cleaning steps:

1. Take out the handpiece, charger, and base on the workbench.

2. Wet the soft cloth completely with distilled water or deionized water, and then wipe all the surfaces of the components such as the handpiece, charger, base, etc. until the surface of the component is not stained.

3. Wipe the surface of the component with a dry soft lint-free cloth until the component is dry.

4. Repeat the above steps at least 3 times.

Note:

a) Use distilled water or deionized water at room temperature for cleaning.

#### 6.3.1.2 Manual disinfection steps:

1. Soak the dry soft cloth with 75% alcohol.

2. Wipe all surfaces of handpiece, charger, base and other components with a wet soft cloth for at least 3 minutes.

3. Wipe the surface of the component with a dry soft lint-free cloth until the component is dry.

Note:

a) The cleaning and disinfection must be performed within 10min before use.

b) The disinfectant used must be used immediately, no foaming is allowed.

c) In addition to 75% alcohol, you can use non-residue disinfectants such as Oxytech from Germany, but you must respect the concentration, temperature and time specified by the disinfectant manufacturer.

d) After cleaning and disinfecting the handpiece, install the disposable sleeve after the machine surface is dry before use, and repeat steps 1, 2 and 3 to clean the disposable insulation sleeve (For detailed installation steps, see section 2.7).

#### 6.3.2 Post-Op processing

After each use, clean and disinfect the handpiece, charger, and base within 30 minutes.

The specific steps are as follows:

Tools: Lint-free soft cloth, tray

1. Remove the contra-angle from the handpiece, place it in a clean tray, and then remove the disposable insulation sleeve from the handpiece.

2. Soak the lint-free soft cloth with distilled water or deionized water, and then wipe all the surfaces of the components such as the handpiece, charger, base, etc. until the surface of the component is not stained.

3. Wet the dry soft cloth with 75% alcohol, and then wipe all surfaces of the handpiece, charger, base and other components for 3 minutes.

4. Put the handpiece, charger, base and other components back into the clean storage area.

Note:

a) The cleaning and disinfection must be performed within 10min before use.

b) The disinfectant used must be used immediately, no foaming is allowed.

c) In addition to 75% alcohol, you can use non-residue disinfectants such as Oxytech from Germany, but you must respect the concentration, temperature and time specified by the disinfectant manufacturer.

### 6.4 The cleaning, disinfection and sterilization of contra-angle, lip hook, file clip, protective silicon cover, touch probe are as follows.

Unless otherwise stated, they will be hereinafter referred to as “products”.

Warning:

The use of strong detergent and disinfectant (alkaline pH > 9 or acid pH < 5) will reduce the life span of products.

In such cases, the manufacturer takes no responsibility.

The products may not be exposed to temperatures above 138°C.

Processing limit

The products have been designed for a large number of sterilization cycles.

The materials used in manufacture were selected accordingly. However with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products.

The maximum number of sterilizations for products is 250 times.

#### 6.4.1 Initial processing

##### 6.4.1.1 Processing principles

It is only possible to carry out effective sterilization after the completion of effective cleaning and disinfection.

Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that the validated parameters are adhered to during every cycle.

Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic, especially with regard to the additional requirements for the inactivation of prions.

##### 6.4.1.2 Post-operative treatment

The post-operative treatment must be carried out immediately, no later than 30 minutes after the completion of the operation. The steps are as follows:

1. Remove the products from the base, and rinse away the dirt on the surface of handpiece with pure water (or distilled water/deionized water);

2. Dry the products with a clean, soft cloth and place them in a clean tray.

Notes:

a) The water used for this must be pure water, distilled water or deionized water.


#### 6.4.2 Preparation before cleaning steps:

Tools: tray, soft brush, clean and dry soft cloth.

1. Remove the shanks/files.

2. Remove the file clip, insulation sleeve, contra-angle and connecting wire from the handpiece in sequence, and then put them into a clean tray;

3. Use a clean soft brush to carefully brush lip hook, file clip, protective silicon cover, touch probe, head and back cover of the contra-angle until the dirt there is no longer any visible dirt on the surface. Then use a soft cloth to dry the products and put them into a clean tray. The cleaning agent can be pure water, distilled water or deionized water.

a) Press the push-button and pull out the shank/file.

b) When removing the protective silicon cover, pull it straight out slowly.

c) When inserting and removing the contra-angle, turn the handpiece power off beforehand.

#### 6.4.3 Cleaning

The cleaning should be performed no later than 24 hours after the operation.

The cleaning can be divided into automated cleaning and manual cleaning.

Automated cleaning is preferred if conditions permit.

##### 6.4.3.1 Automated cleaning

- The cleaner is CE certified in accordance with EN ISO 15883.
- There should be a flushing connector connected to the inner cavity of the product.
- The cleaning procedure is suitable for the product, and the irrigating period is sufficient.

It is recommended to use a washer-disinfector in accordance with EN ISO 15883. For the specific procedure, please refer to the automated disinfection section in the next section "Disinfection".

Notes:

a) The cleaning agent does not have to be pure water. It can be distilled water, deionized water or multi-enzyme.

But please ensure that the selected cleaning agent is compatible with the product.

b) In washing stage, the water temperature should not exceed 45°C, otherwise the protein will coagulate and it would be difficult to remove.

c) After cleaning, the chemical residue should be less than 10 mg/L.

#### 6.4.4 Disinfection

Disinfection must be performed no later than 2 hours after the cleaning phase. Automated disinfection is preferred if conditions permit.

##### 6.4.4.1 Automated disinfection-washer-disinfector

- The washer-disinfector is CE certified in accordance with EN ISO 15883.
- Use high temperature disinfection function. The temperature does not exceed 134°C, and the disinfection at that temperature cannot exceed 20 minutes.
- The disinfection cycle is in accordance with the disinfection cycle as per EN ISO 15883.

Cleaning and disinfecting steps by using a washer-disinfector:

1. Carefully place the product into the disinfection basket. Fixation of product is needed only if there is a risk of the product coming out of the device.

The products are not allowed to come in contact with each other.

2. Use a suitable rinsing adaptor, and connect the internal water lines to the rinsing connection of the washer-disinfector.

3. Start the program.

4. After the program is finished, remove the product from the washer-disinfector, followed by inspection (refer to section "Inspection and Maintenance") and packaging (refer to chapter "Packaging").

Dry the product repeatedly if necessary (refer to section "Drying").

Notes:

a) Before use, you must carefully read the operating instructions provided by the equipment manufacturer to familiarize yourself with the disinfection process and precautions.

b) With this equipment, cleaning, disinfection and drying will be carried out together.

c) Cleaning: (c1) The cleaning procedure should be suitable for the product to be treated. The rinsing period should be sufficient (5 - 10 minutes). Pre-wash for 3 minutes, wash for another 5 minutes, and rinse it twice, with each rinse lasting for 1 minute. (c2) In the washing stage, the water temperature should not exceed 45°C, otherwise the protein will coagulate and it is difficult to remove. (c3) The solution used can be pure water, distilled water, deionized water or multi-enzyme solution, etc., and only freshly prepared solutions can be used. (c4) During the use of cleaner, the concentration and time provided by manufacturer must be observed. The used cleaner is neodisher MediZym (Dr. Weigert).

d) Disinfection: (d1) Direct use after disinfection: temperature 90°C, time 5 min or A0 3000; Sterilize it after disinfection and use: temperature 90°C, time 1 min or A0 600 (d2) For the disinfection here, the temperature is 93°C, the time is 2.5 min, and A0 > 3000

e) Only distilled or deionized water with a low content of microorganisms (<10 cfu/ml) can be used for all rinsing steps. (For example, pure water that is in accordance with the European Pharmacopoeia or the United States Pharmacopoeia).

f) After cleaning, the chemical residue should be less than 10mg / L.

g) The air used for drying must be filtered by HEPA.

h) Regularly repair and inspect the disinfector.

#### 6.4.5 Drying

If your cleaning and disinfection process does not have an automatic drying function, dry it after cleaning and disinfection.

**Methods:**

1. Spread a clean white sheet of paper (white cloth) on the flat table, point the product against the white paper (white cloth), and then dry the product with filtered dry compressed air (maximum pressure 3 bar). The drying cycle is completed when no further liquid is sprayed onto the white paper (white cloth).
2. It can also be dried directly in a medical drying cabinet (or oven). The recommended drying temperature is 80°C ~ 120°C and the time should be 15 ~ 40 minutes.

**Notes:**

- a) The drying of product must be performed in a clean place.
- b) The drying temperature should not exceed 138°C;
- c) The equipment used should be inspected and maintained regularly.

**6.4.6 Inspection and maintenance**
**6.4.6.1 Inspection**


In this chapter, we only check the appearance of the product.

1. Check the product. If there is still visible staining on the product after cleaning/ disinfection, the entire cleaning/ disinfection process must be repeated.
2. Check the product. If it is obviously damaged, smashed, detached, corroded or bent, it must be scrapped and must not be reused.
3. Check the product. If the accessories are found to be damaged, please replace them before use. New accessories for replacement must be cleaned, disinfected and dried.
4. If the service time (number of times) of the product reaches the specified service life (number of times), please replace it in time.

**6.4.6.2 Maintenance**

Sterilizable Oil lubrication shall be applied to dried contra angle.

The nozzle of cleaning lubricant is aligned with the air intake hole at the end of the contra angle to inject oil for 1-2 seconds.

**6.4.7 Packaging**

Install the disinfected and dried product and quickly package it in a medical sterilization bag (FDA cleared wrap or pouch).

**Notes:**

- a) The package used conforms to ISO 11607;
- b) It can withstand high temperature of 138°C and has sufficient steam permeability;
- c) The packaging environment and related tools must be cleaned regularly to ensure cleanliness and prevent the introduction of contaminants;
- d) Avoid contact with parts of different metals when packaging.

**6.4.8 Sterilization**

Use only the following steam sterilization procedures (fractional pre-vacuum procedure\*) for sterilization, and other sterilization procedures are prohibited:

- The steam sterilizer complies with EN13060 or is certified according to EN 285 to comply with EN ISO 17665 or is cleared by FDA.

The validated sterilization cycle is one (1) fraction cycle, three (3) consecutive half cycle and one (1) full cycle.

The parameters are provided in the table below.

Verification of the fundamental suitability of the products for effective steam sterilization was provided by an accredited test laboratory.

**Notes:**

- a) Only products that have been effectively cleaned and disinfected are allowed to be sterilized;
- b) Before using the sterilizer for sterilization, read the Instruction manual provided by the equipment manufacturer and follow the instructions.
- c) Do not use hot air sterilization and radiation sterilization as this may result in damage to the product;
- d) Please use the recommended sterilization procedures for sterilization. It is not recommended to sterilize with other sterilization procedures such as ethylene oxide, formaldehyde and low temperature plasma sterilization. The manufacturer assumes no responsibility for the procedures that have not been recommended. If you use the sterilization procedures that have not been recommended, please adhere to related effective standards and verify the suitability and effectiveness.

\* Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum. The procedure used is to perform steam sterilization through three pre-vacuums.

Mode	Vacuum	Temperature	Cycle	Exposure Time	Drying Time
PreVac	- 80 kPa	(132-134)°C	Fraction	30 seconds	20 minutes
			Half 1	2 minutes	20 minutes
			Half 2	2 minutes	20 minutes
			Half 3	2 minutes	20 minutes
			Full	4 minutes	20 minutes

**6.4.9 Storage**

1. Store in a clean, dry, ventilated, non-corrosive atmosphere with a relative humidity of 10% to 93 %, an atmospheric pressure of 70kPa to 106kPa, and a temperature of -20°C to +55°C;
2. After sterilization, the product should be packaged in a medical sterilization bag or a clean sealing container, and stored in a special storage cabinet. The storage time should not exceed 7 days. If it is exceeded, it should be reprocessed before use.

**Notes:**

- a) The storage environment should be clean and must be disinfected regularly;
- b) Product storage must be batched and marked and recorded.

#### 6.4.10 Transport

1. Prevent excessive shock and vibration during transport, and handle with care;
2. The products should not be mixed with dangerous goods during transport.
3. Avoid exposure to sun or rain or snow during transport.

## 7 Storage, maintenance and transport

### 7.1 Storage

7.1.1 This equipment should be stored in a room where the relative humidity is 10 % ~ 93 %, atmospheric pressure is 70kPa to 106kPa, and the temperature is -20°C ~ +55°C.

7.1.2 Avoid the storage in excessively hot conditions. High temperature will shorten the life of electronic components, damage battery, reshape or melt some plastic.

7.1.3 Avoid the storage in excessively cold conditions. Otherwise, when the temperature of the equipment increases to a normal level, there will be dew that will possibly damage PCB board.

### 7.2 Maintenance

7.2.1 This device does not include spare accessories for repair, the repair should be carried out by authorized person or authorized after-sales center.

7.2.2 Keep the equipment in a dry storage condition during storage.

7.2.3 Protect the equipment from sudden impacts such as shocks or jolts.

7.2.4 Do not soil the equipment.

7.2.5 Calibration is recommended when using a new/different contra angle or after a prolonged period of operation, as the running properties can change with usage, cleaning and sterilization.

7.2.6 Replace the battery if it seems to be running out of power sooner than expected.

### 7.3 Transport

7.3.1 Protect product from impacts, shaking, jolts and shocks during transport. Handle with care and do not transport the product upside down.

7.3.2 Don't put it together with dangerous goods during transport.

7.3.3 Protect from sunlight, rain and snow. Do not allow products to get wet.

## 8 Environmental protection

Please dispose of the product according to the local laws.

## 9 After sales service

From the date this equipment has been sold, based on the warranty card, we will repair this equipment free of charge if there are quality problems.

Please refer to the warranty card for the warranty period.

## 10 European authorized representative

MedNet EC-Rep GmbH  
Borkstrasse 10 · 48163 Muenster · Germany

**Guilin Woodpecker Medical Instrument Co.,Ltd.**

Information Industrial Park, Guilin National High-Tech

Zone, Guilin, Guangxi, 541004 P.R.China

Sales Dept: +86-773-5873196

E-mail: woodpecker@glwoodpecker.com

Website: <http://www.glwoodpecker.com>

## 11 Icons

	Follow Instructions for Use		Serial number
	Date of manufacture		Manufacturer
	Type B applied part		Classe II equipment
<b>IPX0</b>	Ordinary equipment		Recovery
	Used indoor only		Keep dry

	Handle with care		Appliance compliance WEEE directive
	Humidity limitation		Temperature limitation
	Atmospheric pressure for storage		CE marked product
	Authorised Representative in the EUROPEAN COMMUNITY		Medical Device

## 12 Statement

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or reprocessed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

## 13 EMC-Declaration of conformity

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be effected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

### Technical Description Concerning Electromagnetic Emission

Table 1: Declaration - electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions		
The model EnGO-Motor is intended for use in the electromagnetic environment specified below. The customer or the user of the model EnGO-Motor should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model EnGO-Motor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The model EnGO-Motor is suitable for used in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

### Technical Description Concerning Electromagnetic Immunity

Table 2: Guidance & Declaration - electromagnetic immunity

Guidance & Declaration — electromagnetic immunity			
The model EnGO-Motor is intended for use in the electromagnetic environment specified below. The customer or the user of the model EnGO-Motor should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, ±4, ±8, ±15kV air	±8kV contact ±2, ±4, ±8, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for Input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5, ±1kV line to line ±0.5, ±1, ±2kV line to earth	±0.5, ±1kV line to line ±0.5, ±1, ±2kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95% dip in UT.) for 0.5 cycle <5 % UT (>95% dip in UT.) for 1 cycle 70% UT (30% dip in UT) for 25 cycles <5% UT (>95 % dip in UT) for 250 cycles	<5 % UT (>95% dip in UT.) for 0.5 cycle <5 % UT (>95% dip in UT.) for 1 cycle 70% UT (30% dip in UT) for 25 cycles <5% UT (>95 % dip in UT) for 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the models EnGO-Motor requires continued operation during power mains interruptions, it is recommended that the models EnGO-Motor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Table 3: Guidance & Declaration - electromagnetic immunity concerning Conducted RF & Radiated RF

Guidance & Declaration - Electromagnetic immunity			
The model EnGO-Motor is intended for use in the electromagnetic environment specified below. The customer or the user of the models EnGO-Motor should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM frequency band 3 V/m 80 MHz to 2.7 GHz	3V 6V 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the models EnGO-Motor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \times P^{1/2}$ $d = 2 \times P^{1/2}$ $d = 1.2 \times P^{1/2}$ 80 MHz to 800 MHz $d = 2.3 \times P^{1/2}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz end 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model EnGO-Motor is used exceeds the applicable RF compliance level above, the model EnGO-Motor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model EnGO-Motor. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Brasseler®, Komet®, Art2®, CeraBur®, CeraCut®, CeraDrill®, CeraPost®, DC1®, DIA0®, FastFile®, F360®, F6 SkyTaper®, H4MC®, KometBioSeal®, OccluShaper®, OptiPost®, PolyBur®, PrepMarker®, Procodile®, Procodile Q® und SHAX® sind eingetragene Marken der Gebr. Brasseler GmbH & Co. KG.

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## Komet Austria Handelsagentur GmbH