# Let's Go – with EnGO! Step by Step

English, German, French,  
Italian, Spanish

Let's Go – with **EnGO!**  
Step by Step

# Table of contents

- 1 Unboxing
- 2 Basic settings
- 3 Selection of the program
- 4 Modifications to the programs
- 5 Length measurement
- 6 Insertion and removal of the file
- 7 Reprocessing, maintenance and lubrication

# 1 Unboxing

## Components included with the device

1

2

3

4

5

6

7

8

1) EnGO motor + EnGO charging station

2) EnGO AC/DC adapter EU

3) Apex measuring cable

4) EnGO lip clips + measuring rods

The measuring rods can optionally be used with the file clamp during length measurement and facilitate handling, for ex. in case of molars.

5) EnGo file clips

6) EnGo silicone sleeves

For insulating the contra-angle when measuring length.

7) EnGo disposable insulating sheaths

Before each use of the handpiece and after the handpiece has been cleaned and disinfected, fit a disposable insulating sleeve..

8) EnGo nozzle + rubber seal motor

Nozzle for lubricating the contra-angle with oil.  
Replacement rubber seal to protect the motor against the ingress of liquids (e.g. after disinfection with alcohol-soaked cloths).

➤ The shipment also includes instructions for use & Apex Locator trouble shooting sheet

# 2 Basic settings


Detailed information on the EnGo (switching on and off)

Switching on:  
Press the main button to switch on the EnGo.

Switching off:  
Press and hold the P-button while also pressing and holding the main button at the same time.

Variation 2

# 2 Basic settings

## Administration of the basic settings

### Step by Step:

1. The EnGO is switched off.
2. Press and hold the "P" button and while pressing the main button at the same time.
3. The software version appears.
4. The basic settings can now be adjusted.
5. Press the settings buttons "+" / "-" to select from the setting options
6. Press the settings button "P", to confirm your selection and to move on to the next setting
7. Press the main button to confirm all settings.

 The screen switches to standby mode and a program can be selected.

### Variation 1

After how much time should the EnGO switch off automatically?

Auto Power OFF  
5 min

Options to choose from:

- 3 to 30 minutes

After how much time should the EnGO switch to standby mode?

Auto Standby Scr  
30 sec

Options to choose from :

- 3 to 30 seconds

Which is your dominant hand?

Dominant Hand  
Right

Options to choose from

- Right hand
- Left hand

Would you like to activate the calibration function?

Calibration  
OFF

Options to choose from :

- On
- Off

How loud would you like the signal tones of the EnGo to be?

Beeper Volume  
Vol.3

Options to choose from :

- Vol.0 (mute), Vol.1 (quiet), Vol.2 (loud), Vol.3 (very loud)

Would you like to restore the default setting?

Restore Defaults  
OFF

Options to choose from :

- On (yes)
- Off (no)

# 3

## Selection of the program

### Preinstalled programs


On delivery, the EnGO has 10 memory programs (M0-M9) + 5 preinstalled Komet file systems.

#### 10 Memory programs

- The memory location M0 is occupied with EAL – electronic length measurement
- The memory locations M1-M9 are preset with various setting values which can be adjusted to suit individual requirements. (Explanation to follow)


e. g.: Reciprocating with 30° back and forth movement

e. g.: rotary movement at 300 rpm and 2 Ncm

#### 5 preinstalled Komet file systems

- FQ Opener
- FQ Glider
- FQ04 020 – 035
- FQ04 045 – 055
- FQ06

Additional Komet file systems are available (see next page for further details)

# 3

## Selection of the program


Switching between programs + additional Komet systems

### Step by Step:

1. The EnGO is in standby mode.
2. Press the settings buttons "+" / "-" to switch between the preset memory programs.

### Selection of additional Komet file systems:

1. The EnGO is in standby mode. The selected program is irrelevant for this.
2. Press and hold the "P" button for 3 seconds.
3. Further preinstalled Komet systems appear in the list view.
4. You can switch between programs by pressing the settings buttons "+" / "-".
5. Once you have chosen a system, press the main button to confirm your selection.


# 4 Modifications to the programs

Reading and understanding the screen when in standby mode

Selected program:

- Either one of the memory slots M0-M9
- Or one of the preset Komet file systems

Speed setting:  
Rotational speed of the files

Torque setting:  
Torque.

Battery status:

The EnGO battery easily lasts the entire day in the practice, even during the most difficult endo cases.

Operating mode:

- CW: Clockwise rotation
- CCW: Anticlockwise rotation
- RECIP: Reciprocating mode
- ATR: Auto-Torque-Reverse (when a preset torque limit is reached, the motor switches to anticlockwise rotation to release the blocked instrument and reduce the risk of fracture)
- EAL: Length measurement of a root canal

# 4 Modifications to the programs

## Modification of preinstalled memory programs:

### Step by Step:

1. EnGO is in standby mode and the program to be modified is displayed.
2. Press the "P" button.
3. The different modification options can not be changed:

Operation Mode	Speed	Torque Limit	Forward Angle	Reverse Angle	Apical Action	Auto Start	Auto Stop	Flash Bar Position	Apical Slow Down
CW	250 rpm	2.0 Nm	30°	30°	OFF	OFF	OFF		OFF

↳ The modification options are explained on the next few pages. The forward and reverse angle cannot be adjusted in reciprocating mode.

4. You can change the modification buttons with the setting buttons "+" / "-".
5. Press the "P" button to confirm your choice and to automatically move on to the next modification option.
6. Press the main button to confirm all modifications and to apply them to the program.

# 4 Modifications to the programs

Explanation of the modification options

Operation Mode

**CW**

Operation modes

CW: Clockwise rotation of the file (forward rotation)

CCW: Anticlockwise rotation of the file (reverse rotation)

RECIP: Reciprocating mode

ATR: Adaptive auto-torque-reverse

EAL: Electronic length measurement of a root canal (Electronic Apex-Locator)

# 4 Modifications to the programs

Explanation of the modification options

Speed

250 rpm

Rotational speed of the file

In CW mode: Adjustable from 100 rpm – 2,500 rpm

In CCW mode: Adjustable from 100 rpm – 2,500 rpm

In RECIP mode: Adjustable from 100 rpm – 500 rpm

In ATR mode: Adjustable from 100 rpm – 500 rpm

# 4 Modifications to the programs

## Explanation of the modification options

Torque Limit

2.0 Ncm

Torque

In CW mode : Adjustable from 0.4 Ncm – 5.0 Ncm

In CCW mode: Adjustable from 0.4 Ncm – 5.0 Ncm

In RECIP mode: Adjustable from 2.0 Ncm – 5.0 Ncm

In ATR mode: Adjustable from 0.4 Ncm – 4.0 Ncm

# 4 Modifications to the programs

## Explanation of the modification options

Forward Angle

30°

Reverse Angle

30°

Forward and reverse angle

Forward Angle (clockwise)

ATR mode: Adjustable from 60° - 400°

RECIP mode: Adjustable from 20° - 400°

Reverse Angle (anticlockwise)

ATR mode: The reverse angle cannot be larger than the forward angle.

RECIP mode: Adjustable from 20° - 400°

Individually adjustable in 10° steps

# 4 Modifications to the programs

## Explanation of the modification options

Apical Action

OFF

Apical action

Thanks to the integrated length measurement function, the EnGO recognizes when the specified apical reference point has been reached. It automatically switches to reverse mode, stops or switches itself off as required.

**Off:** The function is switched off. The file rotates as before when the reference point has been reached.

**Stop:** The function is switched off. The file automatically stops rotating when the reference point has been reached. It starts rotating again when it is moved slightly upwards, and the blockage is released..

**Reverse:** The file automatically switches into reverse mode when the reference point has been reached or exceeded. It switches back into the previous direction of rotation if it is moved slightly upwards and the blockage is released.

# 4 Modifications to the programs

## Explanation of the modification options

Auto Start

OFF

Automatic start

The file begins to rotate automatically when it is introduced into the canal and more than 2 lines of the canal length bar are illuminated. This function is only available when the apex cable is plugged in.

Off: The function is switched off. The file does not begin to rotate automatically, and the main button has to be used for starting and stopping, as usual.

On: The function is switched on. The motor starts automatically.

# 4 Modifications to the programs

## Explanation of the modification options

Auto Stop

OFF

Automatic stop

The file stops rotating automatically when it is removed from the canal and less than 2 lines on the canal length bar are illuminated. This function is only available when the apex cable is plugged in.

Off: The function is switched off. The file does not stop rotating automatically, and the main button has to be used for starting and stopping, as usual.

On: The function is switched on. The motor stops automatically.

# 4 Modifications to the programs

## Explanation of the modification options

Apical Slow Down

OFF

Apical slow down


The rotation of the file slows down automatically when the tip of the file approaches the reference point.

Off: The function is switched off.

On: The function is switched on.

# 4 Modifications to the programs

Explanation of the modification options

## Setting of the apical reference point

This is the reference point (small arrow) at which the apical action is triggered off.

The apical reference point can be set to anywhere between 2 and AP (apical foramen). (0.5 indicates that the tip of the file is very close to the physiological apical foramen).

The flashing bar indicates roughly how close the file is to the apex.

Please also consult the chapter "length measurement"

# 4 Modifications to the programs

## Modification options IMPORTANT NOTE

Make sure that all settings are adjusted to suit the file to be used. If not, there is a risk that the file may break.

Help with the correct file settings:

The reverse of our blister packs is provided with a label indicating the setting parameters of the file.

Komet Dental · Gebr. Brasseler GmbH & Co. KG  
Trophagener Weg 25 · 32657 Lemgo · Germany

REF FQ04L25 204 025  
LOT 00433712

FQ  
NiTi Feilen  
NiTi files  
Limes NiTi  
Lime NiTi  
Limas NiTi  
L=25mm, 025/.04

Torque 1,8 Ncm  
max. 500 min<sup>-1</sup>  
opt. 300 min<sup>-1</sup>

UDI HBC  
+E226100167731  
/S00433712  
/16D20221121  
/14D20271031\*

STERILE R  
MD  
2022-11-21  
2027-10-31  
Made in Germany  
VS1

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
# 5 Length measurement

Mechanical preparation, combined with length measurement

## Step by Step:

1. Connect the Apex measuring cable to the EnGO via the USB-C socket (on the back of the EnGO handpiece). Use the silicone cover for the contra-angle handpiece. )
2. A lip hook is connected to the light-colored socket of the apex cable. The black socket is not used.
3. The lip clip is attached to the lip of the patient.
4. The bar showing the canal length is displayed on the screen.
5. Automatic functions such as apical action, auto start etc. can be set according to requirement.

# 5 Length measurement

## Use of the motor as apex locator


### Step by Step:

1. Select in the standby mode of the EAL operation mode.
2. Connect the Apex measuring cable to the EnGO via the USB-C socket (on the back of the EnGo handpiece).
3. A lip clip is connected to the light-colored socket of the apex cable.  
The black socket is connected to a file clamp.

The file clip must hold the file correctly. To ensure this, press the button on the file clip with your thumb in the direction of the arrow. Place the holder on the metal upper part of the file and release the button..


4. The lip clip is attached to the patient's lip.
5. The bar showing the canal length is displayed on the screen.

# 5 Length measurement

## Display of the canal measurement

This display appears when a file is in the canal and the lip clip touches the patient's mouth. The bars in the display show the position of the file tip. In EAL mode, the display is enlarged (see illustration on the right) if the distance to the apex is less than 1.0.

The values displayed do not represent the actual length of the apical foramen. They merely show the progression of the file towards the apex.

The digital number "00" indicates that the file has reached the apical foramen.

The digital numbers -1 and -2 indicate that the file has passed the apical foramen.

Subtract 0.5 -1 mm from the measured file length to calculate the working length. These figures are used to estimate the working length of the canal.

# 6 Insertion and removal of the file


## Insertion and removal

### Step by Step (insertion):

1. The EnGO is switched off.
2. Press and hold the push button of the contra-angle (1) and insert the file (3)
3. Turn the file back and forth until it is in line with the inner locking groove and snaps into place.
4. Release the locking button so that the file can engage in the contra-angle handpiece.

### Step by Step (removal):

1. Press the push button on the contra-angle.
2. Pull out the file.

After inserting the file into the contra-angle handpiece, place your hand over the sliding cover to ensure that the file cannot be removed.

Please only use files with shanks that comply with the applicable ISO standard. (ISO standard:  $\text{Ø}2.334 - 2.350 \text{ mm}$ )

# 7


# Reprocessing, maintenance and lubrication

## Reprocessing / cleaning

↳ Reprocessing / cleaning can take place either manually or mechanically.

**i** A Miele adapter for contra-angles with a silicone adapter with a diameter of 16 mm can be used for mechanical cleaning.

A detailed description of the cleaning, disinfection and sterilization procedures can be found in the instructions for use from page 33 onwards.

# 7

# Reprocessing, maintenance and lubrication

## Preprocessing / cleaning

A detailed description of the cleaning, disinfection and sterilization procedures can be found in the instructions for use from page 33 onwards.

1) Cleaning. → Reprocessing / cleaning can take place either manually or mechanically.

A Miele adapter for contra-angle handpieces with a silicone adapter with a diameter of 16 mm can be used for mechanical cleaning.

2) Disinfection.

3) Drying.

4) Inspection.

5) Lubrication. →


Lubrication can be done with the EnGO nozzle included in the delivery in combination with, for example, KaVo lubricating oil

Click on this link for manufacturer details: <https://youtu.be/BgSYIbp0L1s>

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