EnDrive

Endodontic motor

- Ⓔ User manual
- Ⓓ Gebrauchsanleitung
- Ⓕ Mode d'emploi
- Ⓔ Instrucciones de uso
- Ⓘ Istruzioni per l'uso

Based on Manufacturer Rev3 dated 19/06/2024

EnDrive

- 5 – 39 (EN) User manual
- 41 – 75 (DE) Gebrauchsanleitung
- 77 – 111 (FR) Mode d'emploi
- 113 – 147 (ES) Instrucciones de uso
- 150 – 183 (IT) Istruzioni per l'uso

Brasseler®, Komet®, Art2®, CeraBur®, CeraCut®, CeraDrill®, CeraPost®, DC1®, DIAO®, FastFile®, F360®, F6 SkyTaper®, H4MC®, KometBioSeal®, OccluShaper®, OptiPost®, PolyBur®, PrepMarker®, Procodile®, Procodile Q® et SHAX® are registered trademarks of Gebr. Brasseler GmbH & Co. KG.

Some of the products and designations mentioned in the text are trademarked, patented or copyrighted. The absence of a special reference or the sign ® should not be interpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or reprocessed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products are subject to alterations.
Printing errors excepted.

As at August 2024

EnDrive

Endodontic motor

User manual

EnDrive

- 6 Intended use
- 7 Labels
- 8 Symbols
- 9 Warnings, precautions,
contra-indications, adverse reactions
- 14 Step by Step Instructions
- 30 Cleaning, Disinfection
and Sterilization
- 32 Maintenance
- 33 Technical Characteristics
- 36 Product Disposal
- 37 Troubleshooting
- 38 Warranty and service center
- 39 EMC Tables

Intended use

EnDrive is a battery-operated electronic medical device used for the treatment of root pulp infection, which allows root canal preparation by means of continuously or alternately rotating endodontic instruments with torque control, or retreatment of an already treated tooth, with integrated apex locator.

There are no restrictions on the type of patients for use of the device, except patients with active implantable devices such as pacemakers, defibrillators, cochlear implants, on patients allergic to metals, on pregnant or lactating women, or on children. (see contraindications).

The use of this product is intended only for dentists in dental offices, clinics and hospitals. Use of this product by unqualified personnel may result in harm to the patient, operator, or environment.

To ensure the safety, intended functionality and proper maintenance of the **EnDrive**, please follow the instructions given in this user manual. Nextmed accepts no liability if **EnDrive** is used, stored or handled in a manner contrary to the directions in this user manual.

In order to fulfil its intended purpose EnDrive must be used with its file clamp, lip hook and contra-angle handpiece accessories, which can be reused having first reprocessed them as per Chapter 5.

Labels

Main Unit Label

Handpiece Label

Packaging Label

Power Supply Label

Lip Clip Label

File Clamp Label

Contra-angle Label

Symbols

	Medical device		Class II equipment
	Consult instructions for use		BF applied part
	Manufacturer		Temperature range
	Production date		Humidity range
	Alternating voltage		Continuous voltage
	Model		Serial number
	Non-sterile device		CE mark for the battery charger according to LVD 2014/35/EU and EMC Directive 2014/30/EU
	Electric and electronic equipment disposal		Medical device compliance to Regulation 2017/745

Warnings, precautions, contra-indications, adverse reactions

Warnings

Dentist

- Read the manual before using **EnDrive** on a patient. Use **EnDrive** only in accordance with this manual to ensure safe and effective use.
- Ensure to have a suitable working distance to view the display from the user workstation.
- **EnDrive** is intended for endodontic treatment.
- It must only be used by dentists and assistants for cleaning and preparation, and only in accordance with this manual.
- Modifications to the device are not permitted.
- Only use the detachable parts supplied. In particular, only use the supplied handpiece cable and do not connect any other cables to the USB socket on the handpiece.
- Only connect the supplied patient cable and handpiece cable to the main unit.
- Hearing impaired users should familiarise themselves with the various signals before using the device for the first time.

Environmental conditions

- Do not use **EnDrive** in the presence of free oxygen, anaesthetics or flammable substances, or in a potentially explosive environment.
- **EnDrive** must be used and stored in a safe, dry, clean environment.
- Place **EnDrive** on a flat, stable surface. Use **EnDrive** only under the temperature and humidity conditions specified in the specifications in this manual.
- **EnDrive** should not be placed on damp surfaces, wet environments, or near liquids.
- Do not expose **EnDrive** to direct or indirect heat sources.
- Avoid shocks and falls. After a fall, check

that the **EnDrive** components are functioning properly.

- **EnDrive** requires precautions regarding electromagnetic compatibility and must be installed and operated in accordance with the EMC Directive. In particular, **EnDrive** must not be used in the vicinity of fluorescent lamps, radio transmitters, remote controls or RF communication devices even if they comply with the requirements of CISPR11.
- Ensure that the patient cable or handpiece cable does not restrict the operator's mobility during the procedure.
- Do not knot the patient cable or handpiece cable.
- Do not pull on the cables but grasp the connectors to disconnect them. Ensure that the cables are never under tension.
- The use of the **EnDrive** next to another device should be avoided as this can lead to malfunctions. Ensure that EnDrive and the other device, if any, are operating normally.
- Do not expose **EnDrive** to ultrasonic or infrasonic sources.
- Avoid direct light on the screen to ensure full readability.

During treatment

- The user must use their own clinical judgement for appropriate use.
- Ensure that the file selected on the screen is the one actually in use on the contra-angle handpiece. It is recommended to adhere

to the torque and speed settings shown on the screen by default, as recommended by the file manufacturer.

- Do not use reciprocating files with continuously rotating file settings, and vice versa.
- Switch the main unit off via the on/off button in case it does not respond to commands.
- Lubricate the contra-angle handpiece after each sterilisation cycle.
- Do not lubricate the handpiece, as this will damage the motor.
- Make sure to calibrate the contra-angle handpiece after each sterilisation and lubrication cycle using the appropriate menu and not during patient treatment. The calibration process rotates the motor over the entire speed range, so make sure the file is not mounted on the contra-angle during calibration.
- If stopping the motor via the button on the handpiece does not respond, press the symbol of the instrument in use on the main unit.
- For the apex locator function, ensure that the contact point between the endodontic instrument and forceps is kept dry and clean to avoid inaccurate readings.
- The apex locator function during motor treatment only works with the supplied contra-angle handpiece.
- If the file clamp is not properly connected to the endodontic instrument, or is damaged, **EnDrive** may generate incorrect information.

- If the rotary endodontic instrument is not properly connected to the contra-angle handpiece, **EnDrive** may generate incorrect information.
- The use of gloves and a rubber dam is mandatory during treatment.
- Place a sterile disposable film on the **EnDrive** screen to reduce cross contamination.
- Do not spill water or other liquids on **EnDrive** parts.
- Take an x-ray and place the stopper shorter than the expected working length before using **EnDrive**.
- Before starting the apex locator readings, dry the access cavity with cotton wool and ensure that there is no excess fluid inside the tooth canal to be treated. Excess liquid can obstruct readings. Consider that a canal does not have to be completely dry, otherwise the device will not provide readings.
- Be sure to avoid contact between the instrument and any metal restorations or gingiva.
- The **EnDrive** scale does not represent a distinct length or distance in mm or other units of measurement. It only indicates an estimate of the proximity of the instrument tip to the reference point for the working length (apical foramen).
- Be sure to always use metal endodontic instruments.
- Only use endodontic instruments that are undamaged and in good condition, within the working life of the instrument. Refer to the instructions provided by the instrument manufacturer.
- The following patient-related factors may prevent accurate readings:
 - Obstructed canals.
 - Canals with very wide apexes.
 - Root fractures or perforations.
 - Metal crowns or bridges.
 - Restorations and metal fillers.
 - Endodontic instruments that are much smaller than the size of the canal.
 - Endodontic instrument in contact with the gingiva.

Batteries

- Only use the charger supplied. The use of non-original parts can result in safety hazards.
- Check the battery for leaks regularly.
- If fluid leaks from the main unit or the handpiece, the battery may be damaged. Stop using the device immediately. If the leakage occurs during charging, stop charging. Avoid skin contact with the leaked fluid. Contact a service centre for battery replacement.
- Do not recharge the battery if the battery is suspected to be damaged.
- If the handpiece battery does not fully recharge in 3 hours, consider replacing it.
- Ensure that **EnDrive** has sufficient battery on both main unit and handpiece before starting treatment. With a low charge, **EnDrive** should not be used.

Sterilizable accessories

- Lip clip, file clamp and contra-angle hand-piece are supplied non-sterile and must be sterilised before first use and between patients.

Repairs and defects

- Check that **EnDrive** is working properly if it has been subject to shocks or falls. Do not use it if there is any malfunction, breakage or apparent cracks. The safety and effectiveness of the device may be seriously compromised. In the event of a severe shock, **EnDrive** must be inspected by authorised service.
- **EnDrive** contains no user-replaceable parts. Repairs must be carried out by authorised service. Any unauthorised repair, alteration, or modification to **EnDrive** is not permitted and Nextmed is not responsible for such modifications.
- It is mandatory to sterilise sterilisable parts before sending them to a service centre.
- Protect **EnDrive** with appropriate packaging before returning it for service.
- Nextmed can provide the information necessary to repair the device to authorised service centres.

Precautions

- Nextmed disclaims all liability in the event that **EnDrive** is used outside of the indications in this user manual.
- **EnDrive** may not be able to function correctly under all conditions, therefore it is

mandatory to take a preoperative radiograph before using **EnDrive** and compare it with the information provided by the device. In the event of a clinically relevant discrepancy, it is recommended to take a second radiograph with the endodontic instrument positioned at the working length provided by the device.

Adverse reactions

- If **EnDrive** provides incorrect readings and there is no radiographic data for comparison, the following adverse reaction may occur: inadequate root canal treatment due to over- or under-instrumentation.
- If the file is used contrary to the manufacturer's instructions, if it is not new, if the file in use does not match the one selected on the screen, the following adverse reaction may occur: inadequate root canal treatment due to file breakage within the canal.

Contra-indications

- Do not use the **EnDrive** apex locator function on patients with active implantable devices such as pacemakers, defibrillators, cochlear implants.
- Use of **EnDrive** is not recommended on patients with metal allergies.
- Do not use **EnDrive** for implantology or any dental procedure other than endodontics.
- The device is not recommended for use on people who are colour blind or have hearing difficulties unless these difficulties are balanced.

- Safety and efficacy have not been established in pregnant or lactating women or on children.

In case of incidents and/or adverse events derived from the use of the device, contact immediately Nextmed at: +39 0573 535801, info@next-med.net, nextmed@lamiaptec.it.

Step by Step Instructions



4.1 Package Contents.


Carefully remove the device and accessories from the packaging and place them on a flat surface.

Check that **EnDrive** comes with the components listed below and that they are in good condition.

N°	Description	Quantity
1	Main Unit	1
2	Handpiece	1
3	Contra-angle	1
4	Battery Charger	1
5	EU Plug	1
6	UK Plug	1
7	US Plug	1
8	AU Plug	1
9	Patient Cable	1
10	Handpiece Cable	1
11	File Clamp	1
12	Lip Clip	1
13	Instruction for Use	1
14	Bumpers	7

The device can be used both in landscape orientation, with the on/off button on the left side, and in portrait orientation, with the on/off button on the top side.

4.2 Turning on and preparing the system


Before the first use, apply the bumpers to the bottom of the device in the desired orientation. In particular, apply at least 4 bumpers near the corners on the bottom surface in case of landscape orientation and 3 bumpers near the corners on the bottom surface in case of portrait orientation; check for stability. Afterwards recharge the batteries of the main unit and the handpiece.

To recharge the batteries:

- Connect the supplied battery charger to the main unit through the connector provided. Connect the charger to the mains making sure that it conforms to the 100 - 240 V, 50 - 60 Hz range. Disconnect the charger from the mains to isolate the system from the power supply.
- Connect the handpiece cable to the main unit and the handpiece.

By connecting the charger to the main unit, the main unit is turned on.

When the power is turned on, the welcome screen appears, followed by the operating screen:


The charging status icons (see section 4.8) flash during the charging process. When charging is complete, the icons stop flashing and the system is ready for use. Check that the icons indicate the charged battery status.

The external charger and handpiece cable can now be disconnected from the main unit and handpiece respectively, and the system is ready to work.

With the main unit off and charger disconnected, turn it on by holding down the button on the top or left side. The same welcome screens will appear.

With handpiece off and handpiece cable disconnected, turn it on by pressing and holding the on/off button for approximately 1 second:

The handpiece automatically connects to the main unit via Bluetooth. The following icons are selectable on the operating screen:

	Apex locator function
	Motor treatment function
	Starts the motor calibration
	Indication of handpiece battery level
	Indication of main unit battery level
	Indication of the successful connection between main unit and handpiece
	Scrolls pages forward if allowed
	Scrolls pages backward if allowed
	Device settings



4.3 Apex Locator function

Select the correspondent symbol from the menu screen to access the apex locator screen. This mode is used when locating the apex of the root canal without using the motor function:


Attach the lip clip to the patient's lip on the side opposite the tooth being treated and irrigate the canal preferably with NaOCl 5% and wipe off the excess liquid:

Insert an endodontic instrument CE marked of the same size of the apex, similar to the one in the picture below:





into the root canal, and connect the file clamp to the metal part of the instrument, near the handle above the stopper, as in the pictures below:

Move the instrument toward the apical constriction. The device begins to emit sounds of increasing frequency and pitch as the tip of the instrument approaches the apex. The available levels are shown below.


Visual Indication	File progression and sounds	Visual Indication	File progression and sounds
	<p>Coronal and median sections of the canal. Move the endodontic instrument apically. Low-frequency tone.</p>		<p>Middle and apical sections of the canal. Move the endodontic instrument apically. Mid-frequency tone.</p>

Visual Indication	File progression and sounds	Visual Indication	File progression and sounds
	<p>Reference point for determination of working length. High-frequency tone.</p> 		<p>Overinstrumentation. Very high frequency tones.</p> 

In correspondence of the reference point for the length determination it is possible to extract the instrument from the canal, and define the working length for the shaping of the canal, to be performed through the motor function.

4.4 Motor Function

Selecting the correspondent symbol from the menu screen the motor treatment screen appears:

This mode is used to perform root canal shaping using the handpiece and the contra-angle handpiece mounted on it. The handpiece operates without any connection to the main unit via cables. The patient cable can be disconnected from the main unit.


The screen shows the preferred file systems selectable by the operator from the settings menu. A file system is a set of one or more files defined by its manufacturer to shape the root canal. Files within a given file system can be either continuously rotating or reciprocating.

EnDrive has in memory an instrument library containing most of the existing file systems on the market; each file has defined, within its settings the mode of rotation, either continuous or reciprocating.



Insert a CE-marked rotary endodontic instrument with stem diameter 2.334 to 2.35 mm, according to ISO 1797:2017, of the type of this one in the figure below:

on the **contra-angle** in the following way:

Clicking on the name of the desired file system, the operating screen will pop up:

Make sure the handpiece is connected to the main unit via Bluetooth.

Select the endodontic file in use on the screen making sure that manufacturer, conicity and size match.

The file can be selected from those available within the chosen file system using the buttons:

The file in use appears both indicatively as an image:

Torque and speed settings appear in the icon to the side:

The values can be changed via the icon (see settings):

By pressing the on/off button on the handpiece, the motor rotates with the settings of reciprocating motion or continuous rotation, and with the torque and speed values shown on the display. The LED lights up blue when the motor rotates. Through the settings (see section 4.6 Settings menu), two operating modes can be selected for continuous rotation and two for reciprocating function.

In general, the motor can be stopped by pressing on the file icon in use, in case the button on the handpiece does not respond.

Continuous Rotation

- **Autostop:** The motor rotates until it stops at the set torque, then pressing again the on/off button rotates in the opposite direction to disengage the file until the file no longer encounters resistance; then, it rotates forward again. A beep sounds when 70% of the set torque threshold is reached, and the LED on the handpiece flashes.
- **Autoreverse:** The motor rotates until the set torque is reached, then automatically reverses the direction of rotation to disengage the file until the file no longer encounters resistance; then, it rotates forward again. A beep sounds when 70% of the set torque threshold is reached, and the LED on the handpiece flashes.


Reciprocating Rotation

- **Reciproc stop:** the motor performs reciprocating motion and stops when the torque reaches about 6 Ncm. An acoustic signal is emitted when 70% of the set torque threshold is reached, and the LED on the handpiece changes color.
- **Reciproc reverse:** the motor performs reciprocating motion and when it reaches the set torque it automatically switches to continuous rotation to disengage the file. Pressing the on/off button again returns to reciprocating motion. An acoustic signal is emitted when 70% of the set torque threshold is reached and the LED on the handpiece changes color.


4.5 Motor Function with integrated apex locator

To perform apical detection during treatment with motor function (combined measurement), connect the handpiece cable to the main unit and the handpiece. This mode replaces the connection of the file clamp to the hand file as described in section 4.3. The Lip Clip must be applied to the patient, and the directions for use as described in section 4.3 must be followed. The electrical signal is transmitted through the endodontic instrument connected on the contra-angle handpiece provided. During combined apex location with rotating motor, the LED on the handpiece appears blue in the coronal parts of the canal, green in the apical parts, and red in case of over-instrumentation.

Through the settings, the motor stop function can be selected when the level corresponding to the desired level is reached.


Clicking on the dotted line the level at which the motor stops is selected regardless of the torque setting, among the 3 levels of the apex sections:

4.6 Settings Menu

Files parameters Settings

As described in Section 4.4, click on the file settings icon to change the parameters of continuous or reciprocating rotation.

It is recommended not to alter the default settings in the tool library suggested by the file manufacturers, so to not cause any breakage. However, it is possible to change torque, speed, and rotation angles for each file to give the experienced user more flexibility.

Clicking on + and - the relevant torque and speed settings in the ranges of 0,4 - 5 Ncm and 250-1000 rpm can be changed, respectively, and can be saved clicking the check mark. If a file is a reciproc one, it will only be possible to change forward angle, reverse angle, in the ranges 30°- 360°, and speed limited to 3 steps 250, 300, 350 rpm.

press to save the changed values

press to discard the changes

press to restore the default settings of the file

the warning indicates that the settings of the file are different from default

Files Sequence Settings in the File System


Favorite file systems appear in the motor function screen. It is possible to select the favorite files within the various file systems by selecting from the complete library by clicking on the settings icon next to the file system.

The blue symbol next to the file indicates that it has been selected. Only the selected files appear in the favorite systems, as described in section 4.4.

Device Settings

As described in Section 4.2, click on the settings icon to access the device settings.

File Systems

Click here to select the file systems that will appear in the motor operation screen mentioned in Section 4.4.

The blue symbol next to the file system indicates that it has been selected.

Motor Settings

Click here to select the operating modes of the motor as described in Section 4.4.

- Autostop
- Autoreverse
- Reciproc stop
- Reciproc reverse
- Apex stop

The selection of autostop excludes autoreverse and viceversa, and the selection of reciproc stop excludes the selection of reciproc reverse and viceversa.

Volume

Clicking here the volume can be adjusted in 4 steps: 0 - 25% - 50% - 75% - 100%.

SW Update

By clicking here, the device searches for available Wi-Fi networks. By connecting to the preferred one, it is checked if a software update is available and if so, the user is asked to download and install it automatically.

Only software updates that have no significant impact on this user manual are installed automatically, such as bug fixes and minor changes (only software versions with increment of the second digit of software identification X.Y).

Software versions that have a significant impact on the user manual (increment of the first digit of software identification X.Y) can only be done by the service center that provides the updated version of this manual.

Keep the charger connected to the main unit and the handpiece cable connected to the main

unit and handpiece while this operation is performed; without the charger connected, the update is not installed. During this operation, treatment on the patient is not possible.

Check periodically for the availability of software updates.

Default Settings

Clicking here all factory values of the device are restored. Any custom torque, speed and angle settings for each file are deleted. Motor settings are autoreverse for continuous rotation, reciprocate reverse for reciprocating motion. The volume is set at 50%.

4.7 Automatic power off

EnDrive main unit will turn off automatically after 5 minutes from the time of inactivity of the handpiece and inactivity of the user on the main unit itself, unless the main unit is connected to the charger. The handpiece will turn off 5 minutes after the main unit is turned off.

4.8 Battery charging and good practice

The main unit of **EnDrive** is powered by a battery pack made of 10 Nickel Metal hydride cells in series with nominal 12 V voltage.

To recharge the main unit, connect the supplied battery charger to the rear connector. Connect the charger to the mains, making sure it conforms to the 100-240 V, 50-60 Hz range.

The **EnDrive** handpiece is powered by a Li-Ion cell with a nominal voltage of 3,6 V.

To recharge the handpiece, connect it to the main unit via the supplied handpiece cable. The handpiece is recharged only when the motor is not running.

The battery status of both the main unit and the handpiece are indicated on the operation screen. For both parts, the symbology is the same below.

	Main unit battery status.
	Handpiece battery status.
	Battery fully charged.
	Battery partially discharged.
	Battery low, recharge before to start any treatment
	Blinking, battery is charging.

Cleaning, Disinfection and Sterilization

The main unit, handpiece, patient cable, hand-piece cable and charger are not sterilizable. Clean these parts with a commonly used mild detergent.

The File clamp, Lip clip and Contra-angle hand-piece are sterilizable.

Pretreatment: dental pulp and dentin residues should be removed from the accessories immediately; if necessary, use a toothbrush to clean the accessories (within 2 hours of use) with an enzymatic detergent such as CIDEZYME, then rinse. Do not allow residues to dry. Use only clean, soft brushes, a clean, soft cloth or wipe used only for this purpose to manually remove contamination and residue. Do not use metal brushes or straws. Check that no visible contamination or residue remains and, if necessary, repeat the pre-cleaning procedure.

File clamp and Lip clip: clean the parts with CIDEZYME enzymatic detergent, referring to the detergent instructions for proper use and required application time. A detergent or ultra-

sonic bath can be used. Rinse the parts with preferably sterile deionized water. When rinsing the parts, open and close the clamp 4-5 times to help clean the residue.

After cleaning the parts, place them in disinfection bath using Cidex OPA disinfectant for the time prescribed in the manufacturer's instructions for use, then rinse with preferably sterile deionized water. Open and close the clamp 4 – 5 times to facilitate cleaning of debris. Dry with filtered totally oil-free compressed air and make sure the parts are well dried before bagging for sterilization.

Contra-angle: manually clean the outer surfaces with the enzymatic detergent CIDEZYME, referring to the detergent's instructions for proper use and required application time, and inside the endodontic file retaining system (e.g., with a syringe) and rinse with preferably sterile deionized water. After cleaning the surfaces, always disinfect them manually with Cidex OPA for the time prescribed in the manufacturer's instructions for use, then rinse thoroughly with preferably sterile deionized water.

Dry with filtered totally oil-free compressed air and make sure the contra-angle is thoroughly dry before bagging it for sterilization.

Do not immerse the contra-angle handpiece completely in detergent or disinfectant baths or ultrasonic baths!

Insert the previously cleaned and disinfected parts into the sterilization pouches and apply the steam sterilization cycle at 134 °C for 3 minutes. The cycle has been validated for

1,5 minutes. Use autoclaves that comply with EN 13060.


Nextmed does not accept responsibility if these instructions are not followed or if an unvalidated process is used.

Use only properly maintained equipment and personal protective equipment (gloves, goggles) according to national laws and regulations.

Maintenance

Routine maintenance of EnDrive consists of lubricating the contra-angle.

Lubricate the contra-angle after each sterilization cycle with Sirona T1 Spray lubricant as shown:

Spray for about 2 seconds and wipe off oil residue. Rotate the motor a few seconds after lubrication.

Take care that there is no oil residue in the interface part with the motor. Perform the calibration of the contra-angle handpiece after lubrication, and not during patient treatment. The calibration process rotates the motor throughout the speed range, so make sure the file is not mounted on the contra-angle during calibration.

Do not lubricate the motor directly as this will damage the internal gears!

Technical Characteristics

Device

Specification	Description
Manufacturer	Nextmed S.r.l.s Via San Donato 345 51100 Pistoia - Italy
Model	EnDrive
Dimensions	Main Unit: Width: 180 mm Height: 110 mm Depth: 125 mm Handpiece: Length: 200 mm Maximum diameter: 32 mm Minimum diameter: 21 mm
Cable length	Patient cable: approx. 1700 mm Handpiece cable: approx. 1700 mm
Weight	Main unit: 600 g Handpiece: 120 g
Material	Plastic: PC-ABS V-0
Main unit power supply	Battery pack made of 10 AA NiMH cells with PTC fuse Pack capacity: 2000 mAh Rated voltage: 12 Vdc

**Device**

Specification	Description
Handpiece power supply	Li-Ion battery with protection circuit Cell capacity: approximately 3000 mAh Rated voltage: 3.6 V
Main unit power supplier	100-240 V, 50-60 Hz, 24 Vdc, 1,25 A, class II
Applied parts	BF for three applied parts of the same functional group (connections to accessories)
Electrical Safety	Internal battery Class II for the battery charger
Torque range (file)	0,4-5 Ncm (+/-20% at 250 rpm and charged battery)
Speed range (file)	250-1000 rpm (+/-10% at 2 Ncm and charged battery)
Operative conditions	+10/+30 °C, RH<80%
Transport and storage	-20/+70 °C, RH<90%
Protection against liquids and solids	IP20
Classification according to MDR Annex VIII	Ila - Rule 9
Safety level in presence of oxygen or flammable substances	Not suitable
Expected life	5 years
Functioning principle	Apex Locator: via sinusoidal waves at different frequencies Motor: hall sensors brushless
Essential performance	No essential performance

Accessories

Specification	Description
Manufacturer	Nextmed S.r.l.s Via San Donato 345 51100 Pistoia - Italia
Model	File Clamp KOMAPX2 Lip Clip KOMAPX3 Contra-angle KOMCAN
Dimensions	File clamp: length 60 mm Lip clip: diameter 2 mm; length 55 mm Contra-angle: maximum 17 mm; length 66 mm
Cable length	File clamp cable: 60 mm
Weight	File clamp: 6 g Lip clip: 2 g Contra-angle: 30 g
Material	File clamp: nylon e glass fibre Lip clip: stainless steel Contra-angle: anodized aluminium
Reduction ratio	6:1 (Contra-angle)
Operative conditions	+10/+30 °C, RH<80%
Transport and storage	-20/+70 °C, RH<90%
Classification according to MDR Annex VIII	Ila – Rule 5
Expected life	File clamp, Lip clip: 500 sterilization cycles Contra-angle: 1.000 sterilization cycles

Product Disposal

Dispose of the product according to national regulations for end-of-life electronic devices, batteries and potentially contaminated parts.
Do not dispose of in ordinary garbage.

Troubleshooting

Problem	What to do
Unreliable apex locator settings	<p>Check that the endodontic instrument has the size of the root canal.</p> <p>Check that the file clamp is clean and appears golden.</p> <p>Check that the patient cable is not damaged.</p> <p>Check that in case of measurement during motor treatment, the contra-angle handpiece is not damaged; try to check the functionality of apex locator stand-alone by use of use of file clamps.</p>
Motor doesn't run	<p>Disconnect the contra-angle and try the rotation of the motor.</p> <p>If it turns, lubricate the Contra-angle and try again.</p>
Main unit does not charge	Check that the charger is connected and there is power.
Handpiece does not charge	Check that the handpiece cable is not damaged.
Main unit does not turn on	<p>Connect the charger to the mains and to the main unit.</p> <p>Wait 10 minutes and try again.</p>
Handpiece does not turn on	<p>Connect the handpiece cable to the handpiece and the main unit.</p> <p>Wait 10 minutes and try again.</p>
Contra-angle does not calibrate	<p>Disconnect the contra-angle and try the calibration. If it works, lubricate the Contra-angle and try again.</p>
Bluetooth connection doesn't work	<p>Turn off the main unit via the power button without charger connected, then let the handpiece to switch off, and then turn them both on again.</p>

Warranty and service center

The warranty period on the entire device and its parts is 18 months from the date of purchase.

Nextmed Technical Support:

Vi.Vi.Med. S.r.l.

Via Tevere 11, 51018 Pieve a Nievole (PT) -

ITALY

<http://vivimed srl.com/>

EMC Tables

Electromagnetic emissions compliance		
Emission test	Conformity	Electromagnetic environment - guidance
RF Emissions Cispr 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions Cispr 11	Class A	The emissions characteristics of this equipment make it suitable for use in professional healthcare facilities and hospitals. If it is used in a residential environment, this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
Harmonic emissions IEC 61000-3-2	Class A Conforms	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Conform	

Conducted RF IEC 61000-4-6	3 Vrms (0.15 MHz - 80 MHz) 6 Vrms in ISM bands 80% AM at 1 kHz	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the product, including cables. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 V/m 710 MHz, 745 MHz, 780 MHz, 5.240 MHz, 5.550 MHz, 5.785 MHz 27 V/m 385 MHz 28 V/m 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1.720 MHz, 1.845 MHz, 1.970 MHz, 2.450 MHz	
<p>NOTE U_i is the mains voltage in AC before applying test level.</p>		

Electromagnetic immunity compliance			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air		Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4	± 1 kV for I/O lines 100 kHz repetition ± 2 kV for power supply lines 100 kHz repetition frequency		The medical device is suitable for use in professional healthcare establishment environments.
Surge IEC 61000-4-5	± 1 kV line-to-line (± 2 kV line-to-ground; case not applicable to electrically insulated Class II ME Equipment)		
Voltage dips IEC 61000-4-11	0% U _i for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U _i for 1 cycle and 70% U _i for 25/30 cycles. Single phase at 0°		The medical device is suitable for use in professional healthcare establishment environments. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or a battery.
Voltage interruptions IEC 61000-4-11	0% U _i for 250/300 cycles		The medical device is suitable for use in professional healthcare establishment environments.
Power frequency magnetic field IEC 61000-4-8	30 A/m (50 Hz or 60 Hz)		

Recommended minimum separation distances between portable and mobile RF communications equipment and the device as per IEC 60601-1-2 Ed. 3			
<p>The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communication equipment.</p>			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of the transmitter (m)		
	From 150 kHz to 80 MHz d = 1.2 √P	From 80 MHz to 800 MHz d = 1.2 √P	From 800 MHz to 2 GHz d = 1.2 √P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitter rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p>			
<p>NOTES</p> <ul style="list-style-type: none"> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. 			

Brasseler®, Komet®, Art2®, CeraBur®, CeraCut®, CeraDrill®, CeraPost®, DC1®, DIA0®, FastFile®, F360®, F6 SkyTaper®, H4MC®, KometBioSeal®, OccluShaper®, OptiPost®, PolyBur®, PrepMarker®, Procodile®, Procodile Q® und SHAX® sind eingetragene Marken von Gebr. Brasseler GmbH & Co. KG.

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



Indicateur optique	Tracé de la lime et tons	Indicateur optique	Tracé de la lime et tons
<p>Three 3D optical indicators showing the coronal and central sections of a root canal. The indicators are shown in a vertical stack, with the top one showing the coronal section and the bottom one showing the central section.</p>	<p>Sections coronaire et centrales du canal. Déplacer l'instrument endodontique apicalement. Sons à basse fréquence</p> <p>A tracing diagram showing a vertical line representing the instrument's position in the canal. The line is shaded to indicate the sound profile. The sound is characterized by low frequency and low amplitude, with a dashed line indicating the level of the instrument.</p>	<p>Two 3D optical indicators showing the middle and apical sections of a root canal. The indicators are shown in a vertical stack, with the top one showing the middle section and the bottom one showing the apical section.</p>	<p>Sections moyenne et apicale du canal. Déplacer l'instrument endodontique apicalement. Sons de moyenne fréquence.</p> <p>A tracing diagram showing a vertical line representing the instrument's position in the canal. The line is shaded to indicate the sound profile. The sound is characterized by medium frequency and medium amplitude, with a dashed line indicating the level of the instrument. The label "Foramen apicale" is visible at the bottom of the diagram.</p>

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Assistenza tecnica di Nextmed. S.r.l.s: