


Surgical bur block

Height 40 mm / 50 mm / 70 mm



So simple, so convenient

- Universal bur block made of stainless steel, particularly suitable for longer surgical instruments
- Saves space
- Tiltable lid to ensure stable position
- 6 blue silicone plugs as universal instrument holders
- Different types of shanks can be combined. The block fits turbine, hand-piece and contra-angle instruments
- Equally suitable for organizing preparation sequences, e.g. rotary instruments to perform a sinus floor elevation
- Bur block can be reprocessed according to all popular methods (disinfection, ultrasound, thermo disinfectant)



9890L5 (without instruments)

**9890L4 / Height 40 mm suitable
for example for bone cutters:**

H254E/LE.314

Trepan burs:

227A/B/C.204

**9890L5 / Height 50 mm suitable
for example for bone cutters:**

H141A/AZ, H162A/AZ, H166A/AZ.104,

Pilot drills for implantology:

210L16/L19/L20.204 und .205

**9890L7 / Height 70 mm suitable
for instruments with shank 105
long hand-piece**