

Compass | Sonic and ultrasonic tips for periodontal treatments

Recommendations – Products and their use in the dental practice



Removal of subgingival concretions

Lateral tooth surfaces:

- Sonic tips: SF4L/R, SF10L/R
- Sonic tips **Quick**: SFQ4L/R, SFQ10L/R
- Ultrasonic tips **EM1**: PL1, PL2
- Ultrasonic tips **SI1**: PE2, PE3, PE5, PE6
- Ultrasonic tips **KA1**: LU64, LU65

Anterior region:

- Sonic tips: SF4, SF10T
- Sonic tips **Quick**: SFQ4
- Ultrasonic tips **EM1**: PL3
- Ultrasonic tips **SI1**: PE1, 4PS
- Ultrasonic tips **KA1**: LU63




Treatment of furcations

- Ultrasonic tips **EM1**: PL4, PL5
- Sonic tips: SF11
- Sonic tips **Quick**: SFQ24L/R, SFQ26, SFQ27



Root planing

- Sonic tips: SF10L/R/T
- Sonic tips **Quick**: SFQ10L/R
- Ultrasonic tips **EM1**: DPL3, HPL3



Sonic connection

REF	Figure
SF4.000	
SF4L.000	
SF4R.000	
SF10L.000	
SF10R.000	
SF10T.000	
SF11.000	

Sonic tips can be used:

in the sonic hand piece Komet SF1LM/S, in the scalers of the co. W&H (series Proxeo® ZA-55/L/LM/M/LS and Proxeo® ST ZE-55RM/BC, series Synea® ZA-55/L/LM/M or series Alegra® ST ZE-55RM/BC), in the SONICflex™ hand piece of the co. KaVo (series 2000N/L/X/LX or series 2003N/L/X/LX), in the SIROAIR L of the co. Sirona®

Sonic Quick Connection

REF	Figure	
SFQ4.000	SFQ24L.000	
SFQ4L.000	SFQ24R.000	
SFQ4R.000	SFQ26.000	
SFQ10L.000	SFQ27.000	
SFQ10R.000		

The sonic tips Quick can be used:

in the sonic hand piece Komet SFQ2008L/LS, in the sonic hand piece SONICflex™ quick of the co. KaVo, series 2008/S/L/LS

Ultrasonic connection EM1

REF	Figure
PL1.EM1	
PL2.EM1	
PL3.EM1	
PL4.EM1	
PL5.EM1	
DPL3.EM1	
HPL3.EM1	

The ultrasonic tips EM1 can be used:

in the Piezon® Master 700, Minimaster®, Minimaster® LED and MiniPiezon® of the co. EMS, in the Titanus® E of the co. TEKNE DENTAL

Ultrasonic Connection SI1

REF	Figure
PE1.SI1	
PE2.SI1	
PE3.SI1	
4PS.SI1	
PE5.SI1	
PE6.SI1	

The ultrasonic tips SI1 can be used:

in the Periosonic for PerioScan, SIROSONIC/L and SIRONSONIC TL of the co. Sirona®, intended for use in the treatment units TENEO and SINIUS

Ultrasonic connection KA1

REF	Figure
LU63.KA1	
LU64.KA1	
LU65.KA1	

The ultrasonic tips KA1 can be used:

in the PIEZOlux™ or SONOsoft/SONOsoft LUX of the co. KaVo